

Sterile Sampling Bags, LDPE

Keep samples safe—easy-seal bags hold solids and liquids

- Quick and easy sealing—flip bag around wire/tape and bend to seal
- No more poking—puncture-proof tabs keep wire intact even after extensive use

Sampling Bags are available with or without a white labeling area.

“Speci-Sponge” Bags are designed to sample surfaces and equipment for *Listeria*, *Salmonella*, and other food-borne pathogens; contain a 1½" x 3" (3.8 x 7.6 cm) sponge that is free of bactericides and noninhibitory; with white labeling area.

Sodium Thiosulfate Bags are EPA-compliant for potable water collection; contain a nontoxic, nonnutritive 10-mg active ingredient tablet of sodium thiosulfate to neutralize the chlorine at the point of collection; with white labeling area.

High-Temperature Bag is for hot fill samples with a temperature up to 215°F (102°C). The FDA-compliant film is a blend of low-density polyethylene and nylon that will hold liquids, solids, sauces, gravies, and other samples having a high temperature. Bag is self standing and does not need a rack or special holder to keep it upright, making sample collection safe and easy. With white labeling area. **Caution:** Do not put bag in boiling water.



Capacity oz (mL)	Size (W x L)	Thickness (mil)	Catalog number	Qty/pk	Price/pk
Sampling bags without white labeling area					
4 (118)	3" x 7¼" (7.6 x 18.4 cm)	2.25	GY-06499-10	500	
18 (532)	4½" x 9" (11.4 x 22.9 cm)	2.5	GY-06499-20	500	
24 (710)	6" x 9" (15.2 x 22.9 cm)	3.0	GY-06499-30	500	
36 (1065)	5" x 15" (12.7 x 38.1 cm)	3.0	GY-06499-40	500	
42 (1242)	6" x 15" (15.2 x 38.1 cm)	3.0	GY-06499-50	500	
69 (2041)	7½" x 15" (19.1 x 38.1 cm)	3.0	GY-06499-54	500	
Standup sampling bags with white labeling area					
4 (118)	3" x 7¼" (7.6 x 18.4 cm)	2.5	GY-06499-16	500	
18 (532)	4½" x 9" (11.4 x 22.9 cm)	3.0	GY-06499-21	500	
24 (710)	6" x 9" (15.2 x 22.9 cm)	3.0	GY-06499-31	500	
27 (798)	5" x 12" (12.7 x 30.5 cm)	3.0	GY-06499-32	250	
36 (1065)	5" x 15" (12.7 x 38.1 cm)	4.0	GY-06499-33	250	
42 (1242)	6" x 15" (15.2 x 38.1 cm)	4.0	GY-06499-34	250	
69 (2041)	7½" x 15" (19.1 x 38.1 cm)	4.0	GY-06499-35	250	
Sampling bags with white labeling area					
2 (58)	3" x 5" (7.6 x 12.7 cm)	2.5	GY-06499-60	500	
4 (118)	3" x 7¼" (7.6 x 18.4 cm)	2.25	GY-06499-70	500	
7 (207)	3¾" x 7" (9.5 x 17.8 cm)	2.5	GY-06499-72	500	
13 (384)	5¼" x 7½" (13.3 x 19.1 cm)	2.5	GY-06499-74	500	
18 (532)	4½" x 9" (11.4 x 22.9 cm)	2.5	GY-06499-80	500	
24 (710)	6" x 9" (15.2 x 22.9 cm)	3.0	GY-06499-85	500	
27 (798)	5" x 12" (12.7 x 30.5 cm)	3.0	GY-06499-87	500	
69 (2041)	7½" x 15" (19.1 x 38.1 cm)	3.0	GY-06499-89	500	
92 (2721)	10" x 15" (25.4 x 38.1 cm)	4.0	GY-06499-94	250	
123 (3637)	10" x 20" (25.4 x 50.8 cm)	4.0	GY-06499-98	250	
184 (5441)	15" x 20" (38.1 x 50.8 cm)	4.0	GY-06499-99	100	
"Speci-Sponge" bags for surface sampling with sponge					
18 (532)	4½" x 9" (11.4 x 22.9 cm)	2.5	GY-06499-13	100	
55 (1627)	7½" x 11¾" (19.1 x 29.8 cm)	4.0	GY-06499-14	100	
18 (532)	4" x 9" (10.2 x 22.9 cm)	2.5	GY-06499-52†	100	
Sodium thiosulfate bags for water sampling with tablet					
3% (100)	3" x 7¼" (7.6 x 18.4 cm)	2.5	GY-06499-90	100	
10% (300)	4½" x 9" (11.4 x 22.9 cm)	2.5	GY-06499-95	100	
High-temperature bag					
18 (532)	4½" x 9" (11.4 x 22.9 cm)	3.0	GY-06498-16	500	

†Model 06499-52 has a sterile glove.

Accessories

GY-66203-10 Carrying rack for 2- and 4-oz bags. Double-grid design keeps bags upright and has 15 compartments. Vinyl coated, autoclavable

GY-66203-00 Carrying rack for 18-oz bags. Double-grid design keeps bags upright. Rack has 20 compartments. Vinyl coated, autoclavable

GY-06498-12 Water scoop, 18" (45.7 cm) L coated handle. Two alligator clips hold a variety of bag sizes

GY-06498-14 Mailing tube for 3%-oz (100 mL) sodium thiosulfate bag. Polyethylene, 5"L x 2" dia (12.7 x 5.1 cm)



Standup bags stay upright



Water scoop holds bag open for easy filling



Mailing tube



Rack

Zip-Lip Bags, LDPE

Won't let in contaminants

- No need to worry about dust or moisture getting to your dry samples
- Single track seal not for use with liquids
- Easily identify samples by using bags with write-on labels
- FDA-compliant



Size (W x L)		Thickness (mil)	Catalog number	Qty/ pk	Price/ pk
in.	cm				
Zip-lip bags without labeling area					
2 x 3	5.1 x 7.6	2	GY-60101-00	1000	
3 x 5	7.6 x 12.7		GY-60101-10		
4 x 4	10.2 x 10.2		GY-60101-20		
4 x 6	10.2 x 15.2	2	GY-60101-30	1000	
5 x 8	12.7 x 20.3		GY-60101-40		
6 x 6	15.2 x 15.2		GY-60101-50		
6 x 9	15.2 x 22.9	4	GY-60101-60	1000	
8 x 10	20.3 x 25.4	2	GY-60101-70		
9 x 12	22.9 x 30.5	3	GY-60101-80		
12 x 15	30.5 x 38.1	2	GY-60101-90		
Zip-lip bags with labeling area					
2 x 3	5.1 x 7.6	2	GY-60102-00	1000	
3 x 5	7.6 x 12.7		GY-60102-10		
4 x 6	10.2 x 15.2		GY-60102-30		
5 x 8	12.7 x 20.3	2	GY-60102-40	1000	
6 x 9	15.2 x 22.9		GY-60102-60		
Assortment pack, zip-lip bags w/o label area; 100 ea. of 60101-10, -40, -50, -60, & -80					
—	—	—	GY-06503-90	500	

Heavy-Duty Bag, PE

No more weak bags

- Conforms to ordinances for disposal of dead animals
- Ideal for animal research labs, and humane societies



Size (W x L)		Thickness (mil)	Catalog number	Qty/ pk	Price/ pk
in.	cm				
36 x 50	91.4 x 127	4	GY-06504-10	25	

Autoclavable Bags, PP

No worrying about contaminants

- Stands up to your demands—sturdy 1.5-mil polypropylene construction
- Autoclavable at 250°F (121°C)



Size (W x L)		Thickness (mil)	Catalog number	Qty/ pk	Price/ pk
in.	cm				
4 x 6	10.2 x 15.2	1.5	GY-06505-01	100	
4 x 10	10.2 x 25.4	1.5	GY-06505-03	100	
6 x 10	15.2 x 25.4	1.5	GY-06505-05	100	
6 x 12	15.2 x 30.5	1.5	GY-06505-07	100	
10 x 16	25.4 x 40.6	1.5	GY-06505-09	100	
12 x 18	30.5 x 45.7	1.5	GY-06505-11	100	
12 x 24	30.5 x 70	1.5	GY-06505-13	100	
16 x 18	40.6 x 45.7	1.5	GY-06505-15	100	
18 x 23	45.7 x 58.2	1.5	GY-06505-17	100	
22 x 30	55.9 x 76.2	1.5	GY-06505-19	100	

[GY-06831-15](#) Plastic twist ties, 4"L (10.2 cm). Pack of 1000

Bags, LDPE

Stands up to your needs with strong bottom seal

- Seal with portable heat sealer or bag ties



Portable heat sealer 03014-00

Size (W x L)		Thickness (mil)	Catalog number	Qty/ pk	Price/ pk
in.	cm				
3 x 5	7.6 x 12.7	2	GY-06505-02	1000	
3 x 8	7.6 x 20.3		GY-06505-04		
4 x 5	10.2 x 12.7		GY-06505-06		
5 x 7	12.7 x 17.8		GY-06505-08		
6 x 9	15.2 x 22.9	2.5	GY-06505-10	500	
6 x 18	15.2 x 45.7	3	GY-06505-12	500	
8 x 10	20.3 x 25.4		GY-06505-14		
10 x 12	25.4 x 30.5		GY-06505-16		
12 x 18	30.5 x 45.7	2	GY-06505-18	500	
18 x 24	45.7 x 61		GY-06505-20		
18 x 30	45.7 x 76.2		GY-06506-22		

[GY-03014-00](#) Portable heat sealer; 110 VAC, 60 Hz

[GY-03014-02](#) Portable heat sealer; 240 VAC, 50/60 Hz

[GY-06831-15](#) Plastic tie twists, 4"L (10.2 cm). Pack of 1000

Amber Bags, PE

Don't let light get to your samples

- Protect photo-sensitive samples during storage and transport
- Choose resealable zip-top or open-end styles



Size (W x L)		Thickness (mil)	Catalog number	Qty/ pk	Price/ pk
in.	cm				
Resealable zip-top bags					
2 x 3	5.1 x 7.6	3	GY-60103-00	100	
3 x 5	7.6 x 12.7		GY-60103-01		
4 x 6	10.2 x 15.2		GY-60103-02		
5 x 8	12.7 x 20.3		GY-60103-03		
8 x 8	20.3 x 20.3		GY-60103-04		
8 x 14	20.3 x 35.6	3	GY-60103-05	100	
9 x 12	22.9 x 30.5		GY-60103-06		
12 x 12	30.5 x 30.5		GY-60103-07		
3 x 9	7.6 x 22.9		GY-60103-08		
Open-end bags					
5 x 7	12.7 x 17.8	2	GY-60103-09	100	
6 x 10	15.2 x 25.4		GY-60103-11		
8 x 14	20.3 x 35.6		GY-60103-13		
22 x 18	55.9 x 45.7		GY-60103-14		

Bags

Biohazard Disposal Bags, HDPE**Hold more—up to 50 gallons**

- Thick, 1.25-mil HDPE bags stand up to your needs
- No misinterpretations—procedures printed in four languages
- Know when it's autoclaved—indicator patch darkens for verification
- Autoclavable at 250°F (121°C)



Capacity	Size (W x L)		Catalog number	Qty/ cs	Price/ cs
	in.	cm			
2 to 4 gal.	14 x 19	35.6 x 48.3	GY-06494-00	200	
6 to 9 gal.	19 x 23	48.3 x 58.4	GY-06494-02		
10 to 12 gal.	25 x 35	63.5 x 88.9	GY-06494-04		
20 to 30 gal.	31 x 38	78.7 x 96.5	GY-06494-06		
40 to 50 gal.	38 x 48	96.5 x 121.9	GY-06494-08	100	

[GY-06831-15](#) Plastic ties, 4"L (10.2 cm). Pack of 1000**Autoclavable Biohazard Disposal Bags, HDPE****Meet or exceed ASTM D 1709-85**


- Thick, 1.25-mil HDPE construction
- No misinterpretations—procedures printed in four languages
- Know when it's autoclaved—indicator patch darkens for verification
- Autoclavable at 250°F (121°C)



Capacity	Size (W x L)		Catalog number	Qty/ pk	Price/ pk
	in.	cm			
½ to 1 gal.	12 x 24	30.5 x 61.0	GY-06497-30	100	
10 to 12 gal.	24 x 30	61.0 x 76.2	GY-06497-40	100	
15 to 20 gal.	24 x 36	61.0 x 91.4	GY-06497-50	100	

[GY-06497-14](#) Optional stand with tray; for bags up to 9 gal.[GY-06497-18](#) Optional stand with tray; for bags 10 to 12 gal.[GY-06831-15](#) Plastic twist ties, 4"L (10.2 cm). Pack of 1000**Sterile Scoop An' Bag, PS/PE****No need to buy a separate scoop and bag**

- Sterilized and bagged together—lot-traceable packages
- Easy close twirl type closure
- FDA-compliant polystyrene scoop



Scoop capacity mL (oz)	Scoop handle	Bag size (W x L)	Catalog number	Qty/cs	Price/cs
60 (2)	4" (10.2 cm) L	4½" x 9" (11.4 x 22.9 cm)	GY-06260-15	50	

Desktop Biohazard Bag, HDPE**Goes where you need it**

- Versatile bags can be used with optional stand or used to pick up individual items
- Know when it's autoclaved—indicator patch darkens for verification
- Autoclavable at 250°F (121°C)



Size (W x L)		Catalog number	Qty/ pk	Price/ pk	Optional stand	
in.	cm				Cat. no.	Price
8½ x 11	21.6 x 27.9	GY-06496-50	100		GY-06496-60	

[GY-06831-15](#) Plastic twist ties, 4"L (10.2 cm). Pack of 1000**Heavy-Duty Biohazard Disposal Bags, PP****Don't worry about tears and leaks**

- Stands up to your demands—extra strength 2-mil thick polypropylene
- No misinterpretations—procedures printed in four languages
- Know when it's autoclaved—indicator patch darkens for verification
- Autoclavable at 300°F (149°C)



Capacity	Size (W x L)		Catalog number	Qty/ pk	Price/ pk
	in.	cm			
2 to 4 gal.	14 x 19	35.6 x 48.3	GY-06495-00	200	
6 to 9 gal.	19 x 23	48.3 x 58.4	GY-06495-02	200	
10 to 12 gal.	25 x 35	63.5 x 88.9	GY-06495-04	200	

[GY-06497-14](#) Optional stand with tray; for bags up to 9 gal.[GY-06497-18](#) Optional stand with tray; for bags 10 to 12 gal.[GY-06831-15](#) Plastic twist ties, 4"L (10.2 cm). Pack of 1000