

Centrifuge Bottles

Get better separation of large samples

Conical Bottles

- Measure more precisely with molded-in graduations
- Include leakproof, PP cap

Flat-Bottom Bottles

- Use model 06107-20 at low speeds, filled completely, and in refrigerated centrifuge to maintain leak-free spinning
- Bottles 06107-20 and 06106-10 require sealing cap assembly 06327-57 for spinning at high speeds
- 1000 mL bottles require sealing cap 06327-69 for leak-free spinning at maximum speed



Capacity (mL)	Material	Dimensions (H x OD)	Max RCF (x g)	Autoclavable	Bottle			Adapter required		
					Catalog number	Qty/pk	Price/pk	Catalog number	Qty/pk	Price/pk
Conical bottles										
175	PPCO	139.2 x 61.5 mm	27,500	Yes	GY-06107-28	4		GY-06106-78	4	
	PC	140.0 x 61.4 mm			GY-06106-08					
Flat-bottom bottles										
250	PPCO	127.7 x 61.8 mm	13,200	Yes	GY-06105-20	4		GY-17404-06†	6	
	HDPE	131.9 x 61.2 mm	8000	No	GY-06107-10	6				
	FEP	128.8 x 60.0 mm	4000	Yes	GY-06107-20	1				
250	PC	127.6 x 61.8 mm	27,500	Yes	GY-06106-10	4		GY-17404-06†	6	
		135.1 x 61.5 mm			GY-06106-12					
		133.4 x 60.7 mm			GY-06106-14					
450	PPCO	160.0 x 69.5 mm	13,700	Yes	GY-06105-55	4		—	—	—
	PC				GY-06106-35					
500	PPCO	169.8 x 73.8 mm	4800	Yes	GY-06105-50	4		—	—	—
	PPCO	170.2 x 69.5 mm	13,700		GY-06106-25					
	PC	169.6 x 69.5 mm	13,700		GY-06106-30					
1000	PPCO	184.5 x 97.5 mm	7100	Yes	GY-06105-60	4		—	—	—
	PC	188.4 x 97.6 mm			GY-06106-40					
	PPCO	179.0 x 97.7 mm			GY-06107-75					
	PC	180.2 x 98.1 mm			GY-06107-80					

†These bottles require the noted adapter to fit the rotor.

[GY-06327-57](#) Sealing cap assembly, 38-mm. Pack of 2

[GY-06327-69](#) Replacement PP sealing cap for 450, 500, and 1000 mL bottles. Pack of 2

BOD Bottles and Bottle Racks

Easily identify samples

- Clear glass BOD bottles come with glass stoppers and a white write-on marking space

Use a bottle rack for storing and transporting your BOD bottles. Racks are corrosion resistant and plastisol coated. For our full range of benchtop dissolved oxygen meters including BOD, see your local representative.

Specifications

Capacity: 300 mL

Description	Number sequence	Qty/ cs	Catalog number	Price/ cs
Penny head stopper BOD bottles, clear	Unnumbered	24	GY-03660-14	
Robotic stopper BOD bottles, clear	Unnumbered	24	GY-03660-15	
Robotic stopper BOD bottles, clear	1–24	24	GY-03660-16	
	25–48		GY-03660-17	
	49–72		GY-03660-18	
	73–96		GY-03660-19	



Accessories

[GY-03660-09](#) Cap and insert for use with 03660-15 BOD bottles. Pack of 50

[GY-03660-12](#) Rubber bottle leash holds stopper to bottle. Pack of 25

[GY-03660-20](#) BOD bottle rack, holds twelve 300-mL bottles

Safety Bottles

Protect yourself and your samples

Amber and Clear Safety Bottles

Plain amber bottles (USP Type III) are available in 250-, 1000-, and 2500-mL sizes. PVC-coated amber bottle 07914-84 is designed for safety—won't shatter like regular bottles. Clear, borosilicate glass bottle 07914-82 features a sidearm for a 0.2- μ m air filter.



Color/type	Capacity (mL)	Screw cap size (mm)	Catalog number	Price/ea
Amber/square	250	33	GY-07914-80	
Clear/with side arm	500		GY-07914-82	
Amber/PVC coated	500		GY-07914-84	
Amber/square	1000	45	GY-07914-86	
Amber/round	2500		GY-07914-88	

Safety Bottles with PTFE-Lined Caps

Safety glass bottles are coated to help you meet OSHA requirements by reducing breakage. The plastisol coating helps protect you from acids, solvents, and chemicals. Bottles include caps with a PTFE lining.



Color	Capacity	Screw cap size (mm)	Catalog number	Qty/cs	Price/cs
Clear	16 oz	28–400	GY-08922-50	24	
	32 oz	33–400	GY-08922-60	12	
Amber	16 oz	28–400	GY-08922-55	24	
	32 oz	33–400	GY-08922-65	12	
	1 gal.	38–430	GY-08922-75	4	

Amber Glass Dropper Bottles

Don't let light harm contents

- Protect light-sensitive materials
- Prevent rolling with flat side
- Meet ASTM E 438, Type II and USP Type III standards
- Includes 1-mL dropper



Capacity mL (oz)	Dropper material	Catalog number	Qty/cs	Price/cs
30 (1)	Glass	GY-34523-38	12	
30 (1)	Plastic	GY-34523-48	72	
60 (2)		GY-34523-50	12	

Graduated Square Bottles, PC

Stop drips as you pour

- Drip-proof spout for better control of liquid
- Know your volume with molded-in graduations
- Autoclavable



Capacity mL (oz)	Cap size	Catalog number	Qty/pk	Price/pk
125 (4)	38-430	GY-06040-50	6	
250 (8)	38-430	GY-06040-55	6	
500 (16)	38-430	GY-06040-60	4	
1000 (32)	38-430	GY-06040-65	4	
2000 (64)	53B	GY-06040-70	1	
Replacement caps, PP				
—	38-430	GY-63000-18	12	
—	53B	GY-63000-20	12	

Dropper Bottles, LDPE

Easy dispensing of solutions

- Keep liquid dispensing controlled with easy dispensing cap
- 15-mL bottle has 2-piece dispensing unit, large sizes have hinge caps



Capacity mL (oz)	Catalog number	Qty/pk	Price/pk
15 (½)	GY-62202-10	12	
30 (1)	GY-62202-20		
60 (2)	GY-62202-30		
125 (4)	GY-62202-40	12	
175 (6)	GY-62202-50		
250 (8)	GY-62202-60		
500 (16)	GY-62202-70		

Environmental Express® Disposable BOD Bottles

Reduce contamination and increase test accuracy

NEW

- No precleaning required!
- Laboratory tested and US EPA approved
- Unbreakable and recyclable
- Packed in plastic-lined boxes to eliminate dust and residue



ENVIRONMENTAL EXPRESS®

Capacity	Catalog number	Qty/pk	Price/pk
300 mL	GY- 53200-02	100	

[GY-53200-12](#) BOD bottle stopper, acrylic. Pack of 25

[GY-53200-14](#) BOD overcap prevents evaporation. Pack of 50

Bottles

1395 Series Media Bottles

Permanent white graduations won't wear off

- Use for mixing or storing of your solutions—full capacity indicated by bead line on 100-mL to 2-L sizes
- Linerless, autoclavable polypropylene screw cap included on each bottle
- Color code your samples with optional color-coded PP caps

Caution: Bottles larger than 2 liters should not be used under vacuum pressure, as breakage may occur. Do not tighten caps immediately after autoclaving.



Capacity mL (oz)	Graduation range (mL)	Thread size	Neck ID	Qty/ cs	Catalog number	Price/ cs
25 (¾)	10 to 25	GL25	14.5 mm	10	GY-34514-20	
50 (1½)	20 to 50	GL32	16.5 mm	10	GY-34514-21	
100 (3)	40 to 80	GL45	29 mm	10	GY-34514-22	
250 (8)	50 to 200	GL45	29 mm	10	GY-34514-23	
500 (16)	100 to 400	GL45	29 mm	10	GY-34514-24	
1000 (32)	100 to 900	GL45	29 mm	10	GY-34514-25	
2000 (64)	400 to 1800	GL45	29 mm	4	GY-34514-26	
5000 (160)	500 to 4500	GL45	29 mm	1	GY-34514-28	
10,000 (320)	2000 to 9000	GL45	29 mm	1	GY-34514-30	

[GY-34514-50](#) Polypropylene cap, orange, GL45. Cs of 20

[GY-34514-53](#) Polypropylene cap, green, GL45. Cs of 20

[GY-34514-54](#) Polypropylene cap, gray membrane, GL45. Cs of 10

[GY-34514-55](#) Polypropylene cap, red high-temp. GL45. Cs of 10

[GY-34514-60](#) Pouring ring, clear, low temperature. Cs of 50

[GY-34514-65](#) Pouring ring, red, high temperature. Cs of 50

51395 Series Media Bottles

Safely store light-sensitive materials

- Store and transport light-sensitive solutions in the range of 3000 to 5000 Å
- Will not drip while pouring
- Don't overfill with beaded fill line
- Autoclavable bottle and PP lid



Capacity mL (oz)	Graduation range (mL)	Graduation interval (mL)	Thread size	Neck ID	Qty/ cs	Catalog number	Price/ cs
100 (3)	40 to 80	20	GL45	29 mm	4	GY-34574-21	
250 (8)	50 to 200	50				GY-34574-22	
500 (16)	100 to 400	100				GY-34574-23	
1000 (32)	100 to 900	100				GY-34574-24	

Media Bottles

Easy, drip-free pouring

- Stop worrying about not knowing volume—bottles have permanent white graduations and marking spot
- Use up to 840°F (450°C)
- Autoclavable bottle and PP lid



Capacity mL (oz)	Graduation range (mL)	Graduation interval (mL)	Thread size	Qty/ cs	Catalog number	Price/ cs
100 (3)	20 to 80	10	GL45	10	GY-34523-00	
250 (8)	50 to 200	25	GL45	10	GY-34523-02	
500 (16)	100 to 400	50	GL45	10	GY-34523-04	
1000 (32)	100 to 900	50	GL45	10	GY-34523-06	
2000 (64)	400 to 1800	100	GL45	4	GY-34523-08	
5000 (160)	500 to 4500	250	GL45	1	GY-34523-10	
10,000 (320)	2000 to 9000	500	GL45	1	GY-34523-12	

[GY-34523-11](#) Replacement cap, PP, blue. Pack of 10

[GY-34523-13](#) Replacement pouring ring, PP, clear. Pack of 10

DON'T RISK IT!

Our chemical compatibility database provides important information on chemical resistance—check it out at ColeParmer.com/ChemChart



Amber Media Bottles

Stop UV-light from affecting samples

- Tinted glass blocks more than 99% of UV light below 400 nm
- Permanent graduations will not wear off
- Dripless pouring with included pouring ring
- Autoclavable lid and bottle



Capacity mL (oz)	Graduation range (mL)	Graduation interval (mL)	Thread size	Qty/ cs	Catalog number	Price/ cs
500 (16.9)	100 to 400	50	GL45	1	GY-34523-20	
1000 (33.8)	100 to 900	50	GL45	1	GY-34523-22	
2000 (67.6)	400 to 1800	100			GY-34523-24	
5000 (169.0)	500 to 4500	250			GY-34523-26	

[GY-34523-11](#) Replacement cap, PP, blue. Pack of 10

[GY-34523-13](#) Replacement pouring ring, PP, clear. Pack of 10

61626 Series Media Bottles

Prevent shattered glass in your facility

- Protective PVC coating for safer handling
- Resists thermal shock; should not be used over direct heat or flame
- Dripless pouring ring included
- Autoclavable bottle and PP cap



Capacity mL (oz)	Graduation range (mL)	Graduation interval (mL)	Thread size	Qty/ cs	Catalog number	Price/ cs
100 (3)	40 to 80	20	GL45	4	GY-34575-21	
250 (8)	50 to 200	50			GY-34575-22	
500 (17)	100 to 400	100			GY-34575-23	
1000 (32)	100 to 900	100	GL45	4	GY-34575-24	
2000 (64)	400 to 1800	200			GY-34575-25	

Narrow-Mouth Bottles, HDPE

Ship samples without worry

- Chemical-resistant material stands up to most corrosives
- Will not leak with PP screw closures



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
4 (1/8)	13 mm	12	GY-06035-02	
8 (1/4)	20 mm		GY-06035-03	
15 (1/2)	20 mm		GY-06035-05	
30 (1)	20 mm	12	GY-06035-10	
60 (2)	20 mm		GY-06035-20	
125 (4)	24 mm		GY-06035-30	
175 (6)	24 mm	12	GY-06035-40	
250 (8)	24 mm		GY-06035-50	
500 (16)	28 mm		GY-06035-60	
500 (16)	38–430	6	GY-06035-80	
1000 (32)	38–430		GY-06035-70	

Cole-Parmer® Class A Glass Media Bottles

Easy pouring with wide-mouth opening

NEW

- Ideal for mixing, storing, and sampling
- Get quick identification—large marking spot and gradations
- Autoclavable PP cap and bottle



Cole-Parmer®

Capacity mL (oz)	Graduation range (mL)	Graduation interval (mL)	Thread size	Qty/ cs	Catalog number	Price/ cs
100 (3)	20 to 80	25	GL45	10	GY-34505-72	
250 (8)	50 to 200	50		10	GY-34505-74	
500 (16)	100 to 400	100		10	GY-34505-76	
1000 (32)	100 to 900	100	GL45	10	GY-34505-78	
2000 (64)	400 to 1800	100		4	GY-34505-80	

Safety-Coated Media Bottles

No more cleaning up shattered glass

- Plastic coating prevents shattering for safer transport
- Easy pouring, wide opening
- Permanent, white gradations and marking spot
- Autoclavable bottle and cap



Capacity mL (oz)	Graduation range (mL)	Graduation interval (mL)	Thread size	Qty/ cs	Catalog number	Price/ cs
100 (3)	20-80	10	GL45	4	GY-34523-15	
250 (8)	50-200	25	GL45	4	GY-34523-17	
500 (16)	100-400	50	GL45	4	GY-34523-19	
1000 (32)	100-900	50	GL45	4	GY-34523-21	
2000 (64)	400-1800	100	GL45	4	GY-34523-23	
5000 (160)	1000-4000	250	GL45	1	GY-34523-25	
10000 (320)	2000-8000	500	GL45	1	GY-34523-27	

Cole-Parmer® Economical Bottle Carrier

Easily and safely transport solutions

- Lightweight, durable carrier transports up to six bottles measuring 3 1/2" diameter or less

Cole-Parmer®



Hole diameter	Qty/pk	Catalog number	Price/pk
3 1/2" (89 mm)	1	GY-06722-04	

Bottles

Cole-Parmer® Narrow-Mouth Bottles,
HDPE, LDPE, PP**Keep samples contained with leakproof design**

- Leakproof PP closure for safer transport
- HDPE bottles are FDA-compliant
- LDPE bottles feature transparent, food-grade design
- PP bottles offer superior chemical resistance and are autoclavable



“I will continue to buy this product.”



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk	Qty/ cs	Catalog number	Price/ cs
HDPE							
30 (1)	22 mm	12	GY-06038-61		72	GY-06038-62	
60 (2)	22 mm	12	GY-06038-63		72	GY-06038-64	
125 (4)	28 mm	12	GY-06038-65		72	GY-06038-66	
250 (8)	28 mm	12	GY-06038-67		72	GY-06038-68	
500 (16)	28 mm	12	GY-06038-69		48	GY-06038-70	
1000 (32)	28 mm	6	GY-06038-71		24	GY-06038-72	
2000 (64)	45 mm	1	GY-06038-73		6	GY-06038-74	
LDPE							
60 (2)	22 mm	12	GY-06038-21		72	GY-06038-22	
125 (4)	28 mm	12	GY-06038-23		72	GY-06038-24	
250 (8)	28 mm	12	GY-06038-25		72	GY-06038-26	
500 (16)	28 mm	12	GY-06038-27		48	GY-06038-28	
1000 (32)	28 mm	6	GY-06038-29		24	GY-06038-30	
PP Autoclavable							
30 (1)	22 mm	12	GY-06038-41		72	GY-06038-42	
60 (2)	22 mm	12	GY-06038-43		72	GY-06038-44	
125 (4)	28 mm	12	GY-06038-45		72	GY-06038-46	
250 (8)	28 mm	12	GY-06038-47		72	GY-06038-48	
500 (16)	28 mm	12	GY-06038-49		48	GY-06038-50	
1000 (32)	28 mm	6	GY-06038-51		24	GY-06038-52	
2000 (64)	45 mm	1	GY-06038-53		6	GY-06038-54	

**1396 Series
Square Bottles****Maximize storage space**

- Square design provides easier storage with no wasted space between bottles
- Drip-free pouring
- Autoclavable bottle and PP cap



Capacity mL (oz)	Thread size	Qty/ cs	Catalog number	Price/ cs
100 (3)	GL32	10	GY-34594-22	
250 (8)	GL45	10	GY-34594-24	
500 (16)	GL45	10	GY-34594-26	
1000 (32)	GL45	10	GY-34594-28	

Narrow-Mouth Bottles, LDPE
Won't break or leak

- See samples more easily with translucent material
- Very low metal content, ideal for trace metal analysis
- Include leakproof, PP caps



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
8 (1/4)	20 mm	12	GY-06033-03	
15 (1/2)			GY-06033-05	
30 (1)			GY-06033-10	
60 (2)			GY-06033-20	
125 (4)	24 mm	12	GY-06033-30	
250 (8)	24 mm		GY-06033-50	
500 (16)	28 mm		GY-06033-60	
500 (16)	38–430		GY-06033-80	
1000 (32)	38–430	6	GY-06033-70	
Replacement caps, PP				
—	24 mm	12	GY-63000-04	
—	28 mm		GY-63000-06	
—	38–430		GY-63000-18	

**Cole-Parmer®
Round Jugs, HDPE****Easier carrying with
molded-in handle**

- Lightweight design won't add unnecessary weight
- Leak-free polypropylene (PP) cap

**NEW**

Capacity L (gal.)	Qty/ pk	Catalog number	Price/ pk
2 (1/2)	24	GY-06089-51	
4 (1)	12	GY-06089-53	

DON'T RISK IT!

Our chemical compatibility database provides important information on chemical resistance—check it out at ColeParmer.com/ChemChart



Narrow-Mouth Bottles, PP

Store strong chemicals without concern

- Bottles have excellent chemical resistance
- Meet requirements of Food Additive Amendments of the Federal Food, Drug, and Cosmetic Act
- Will not leak
- Autoclavable bottle and PP cap



Amber Narrow-Mouth Bottles, HDPE

Protect samples in shipping and storage

- Reduce UV light transmission
- Meet USP requirements for storage of light-sensitive materials
- Leakproof polypropylene (PP) screw closures



Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
4 (1/8)	13 mm	12	GY-06034-02	
8 (1/4)	20 mm		GY-06034-03	
15 (1/2)	20 mm		GY-06034-05	
30 (1)	20 mm		GY-06034-10	
60 (2)	20 mm	12	GY-06034-20	
125 (4)	24 mm		GY-06034-40	
250 (16)	24 mm		GY-06034-50	
500 (32)	28 mm		GY-06034-60	
1000 (64)	38-430	6	GY-06034-70	
Replacement caps, PP				
—	24 mm	12	GY-63000-04	
—	28 mm		GY-63000-06	
—	38-430		GY-63000-18	

Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
4 (1/8)	13 mm	12	GY-06026-02	
8 (1/4)	20 mm		GY-06026-03	
15 (1/2)	20 mm		GY-06026-05	
30 (1)	20 mm		GY-06026-10	
60 (2)	20 mm	1	GY-06026-20	
125 (4)	24 mm		GY-06026-40	
250 (8)	24 mm		GY-06026-50	
500 (16)	28 mm		GY-06026-60	
1000 (32)	38-430	6	GY-06026-70	

Economy Sample Bottles

Identify samples quickly and easily

- Large writing area allows for easy identification
- USP Type I borosilicate glass maintains purity of contents
- Designed to make pouring and cleaning easy
- Graduated to approximate volume
- PTFE-lined screw caps
- Autoclavable



Narrow-Mouth Gray Graduated Bottles, PE

Protect and store photosensitive materials

- Measure more accurately with molded-in graduations
- Easy-grip PP screw cap features tag holes that safely fasten to bottle shoulder holes
- Removable LDPE neck insert provides additional security
- Decrease residue deposits with sloping shoulder design



Capacity mL (oz)	Graduation intervals (mL)	Cap size	Qty/cs	Catalog number	Price/cs
250 (8)	25	33-430	48	GY-08913-25	
500 (16)	50		24	GY-08913-45	
1000 (32)	50		24	GY-08913-55	

Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
50 (1.7)	18.5 mm	50	GY-06035-82	
125 (4)	18.5 mm	50	GY-06035-84	
250 (8)	23 mm	50	GY-06035-86	
500 (16)	23 mm	25	GY-06035-88	
1000 (32)	33 mm	10	GY-06035-90	
2000 (64)	33 mm	2	GY-06035-92	

Bottles

Heavy-Duty Bottles, PP

Withstand repeated full vacuum applications

- Use under repeated full vacuum for 24 hours at 68°F (20°C)
- PP closure with thermoplastic elastomer (TPE) gasket
- Autoclavable



Capacity, mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
250 (8)	53B	6	GY-06257-05	
1000 (32)		6	GY-06257-10	
2000 (64)		2	GY-06257-20	
4000 (128)	83B	1	GY-06257-00	
5000 (160)			GY-06257-30	

Narrow-Mouth Bottles, PFA

Stop worrying about contaminants

- Noncontaminating material has lowest metal ion content of all fluoropolymers
- Ideal for environmental sampling
- Leakproof, heat-resistant design
- Autoclavable



Capacity, mL (oz)	Graduations	Neck size (mm)	Catalog number	Price/ea
100 (3)	50 mL	16.7	GY-06040-30	
250 (8)	50 mL	16.5	GY-06040-32	
500 (16)	100 mL	17	GY-06040-34	
1000 (32)	100 mL	25.3	GY-06040-36	
3000 (96)	500 mL	36.6	GY-06040-40	
5000 (160)	1000 mL	36.4	GY-06040-42	
10,000 (320)	None	36.4	GY-06040-44	

Cole-Parmer® Square One-Quart Tall and Short Bottles, HDPE

Meet FDA-compliance

- Leakproof, PP cap



Description	Cap size	Qty/pk	Catalog number	Price/pk
Square tall	38 mm	36	GY-62270-03	
Square short, with handle		54	GY-62270-13	

**Heavy-Duty Bottles, HDPE**

Strong, leakproof storage

- PP closure will not leak
- Stand up to general storage or shipping



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/pk
1000 (32)	53B	6	GY-06025-50	
2000 (64)	53B	2	GY-06025-60	
4000 (128)	83B	1	GY-06025-70	
Replacement caps, PP				
—	53B	12	GY-63000-20	
—	83B	2	GY-06258-30	

DON'T RISK IT!

Our chemical compatibility database provides important information on chemical resistance—check it out at ColeParmer.com/ChemChart

**Filling/Venting Caps, PP**

Quick and easy aseptic liquid transfer

Filling/Venting Caps have thermoplastic elastomer (TPE) gaskets to allow aseptic liquid transfer. Include 6" (15 cm) and 30" (76 cm) lengths of platinum-cured silicone tubing. Autoclavable.

Fitting Kits include two PP fittings, two acetal nuts, two silicone gaskets, thermoplastic elastomer (TPE) port cap, instructions, and template. Autoclavable.



Tubing ID	Ports	Cap size	Qty/pk	Catalog number	Price/pk
Filling/venting caps					
1/2"	2 tubing, 1 venting	83B	1	GY-06258-00	
1/4"		53B	1	GY-06258-10	
1/4"		83B	1	GY-06258-20	
Barbed bulkhead fitting kits					
1/4"	—	—	2	GY-06259-10	
1/2"	—	—	2	GY-06259-00	

Narrow-Mouth PassPort™ IP2 Bottles, HDPE

Ship hazardous
samples reliably

- Safe for small quantities of corrosives or cleaning supplies
- Leakproof, PP screw cap
- Volume, registration notch, material, and IP2 markings molded into bottom



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
60 (2)	20 mm	12	GY-06253-10	
125 (4)	24 mm	12	GY-06253-20	
250 (8)	24 mm	12	GY-06253-30	
500 (16)	28 mm	12	GY-06253-40	
1000 (32)	38-430	6	GY-06253-50	
4000 (128)	38-430	1	GY-06253-60	

Large Narrow-Mouth Bottles, PP

Store large samples
of chemicals safely

- Leakproof, polypropylene caps
- Built-in shoulder loop for easier handling on 4- and 8-liter sizes
- Autoclavable



Capacity L (gal.)	Cap size	Qty/ pk	Catalog number	Price/ pk
Clear bottles, PP				
2 (½)	38-430	1	GY-06028-10	
4 (1)	38-430		GY-06028-20	
8 (2)	53B		GY-06028-30	
Amber bottles, PP				
2 (½)	38-430	1	GY-06027-10	
4 (1)	38-430		GY-06027-20	
8 (2)	53B		GY-06027-30	

Jugs with Handle, LDPE and PP

Easy pour design

- Carry samples easily with molded-in handle
- PP versions are autoclavable
- Leakproof



Capacity L (gal.)	Cap size	Qty/ pk	Catalog number	Price/ pk
Jugs, LDPE				
4 (1) 8 (2)	38-430 53B	1	GY-06062-10 GY-06062-20	
Jugs, PP, Autoclavable				
4 (1) 8 (2)	38-430 53B	1	GY-06064-10 GY-06064-20	
Replacement caps, PP				
— —	38-430 53B	12	GY-63000-18 GY-63000-20	

Sanitary Bottles, 316L SS (Light Gauge)

Keep air and
liquids out

- Constructed of 316L stainless steel with #9, 24 Ra finish
- Included sanitary Tri-Clamp® flange opening, end cap platinum-cured silicone gasket, and heavy-duty clamp



Capacity L (gal.)	Tri-Clamp opening	Neck ID	Bottom OD	Qty/ pk	Catalog number	Price/ pk
1 (¼)	2" 3"	1⅞" 3"	4" (10.2 cm)	1	GY-07506-01 GY-07506-31	
2 (½)	2" 4"	1⅞" 3⅞"	4¾" (12.1 cm)	1	GY-07506-02 GY-07506-32	
5 (1⅓)	2" 4"	1⅞" 3⅞"	7⅞" (18.1 cm) 6⅝" (16.2 cm)	1	GY-07506-05 GY-07506-35	
10 (2½)	2" 4"	1⅞" 3⅞"	8¾" (22.2 cm) 8⅝" (21.3 cm)	1	GY-07506-10 GY-07506-80	
20 (5¼)	2" 4"	1⅞" 3⅞"	10⅝" (27.0 cm)	1	GY-07506-20 GY-07506-90	

Cole-Parmer® Precleaned EPA Sample Containers**No contaminants—precleaned and ready to use**

- Include certificate of conformance to EPA specifications
- Identify samples easier with included labels (unattached)

A-D. Glass Bottles, Jars, and Jugs are processed and certified to meet the requirements of the current US EPA “Specifications and Guidance for Contaminant-Free Sample Containers” document for semi-volatile organics, pesticides, PCBs, metals, cyanide, fluoride, nitrate and nitrite analyses. Include PTFE-lined caps.

E. Glass Septa Jars are processed and certified to meet the requirements of the current US EPA “Specifications and Guidance for Contaminant-Free Sample Containers” document for purgeable volatile organic analysis. Include PTFE-lined silicone septa bonded to open-top caps.




Capacity mL (oz)	Neck size (mm)	Qty/ cs	Catalog number	Price/ cs
A. Clear wide-mouth glass jars				
60 (2)	53 mm	24	GY-99540-40	
125 (4)	58 mm	24	GY-99540-41	
250 (8)	70 mm	12	GY-99540-42	
500 (16)	89 mm	12	GY-99540-43	
950 (32)	89 mm	12	GY-99540-44	
B. Amber narrow-mouth glass bottles				
125 (4)	22 mm	24	GY-99540-30	
250 (8)	24-414	12	GY-99540-31	
500 (16)	28 mm	12	GY-99540-32	
1000 (34)	33-430	12	GY-99540-33	
C. Amber wide-mouth glass bottles				
125 (4)	38 mm	12	GY-99540-51	
250 (8)	45 mm	12	GY-99540-53	
500 (16)	53 mm	12	GY-99540-55	
950 (32)	53 mm	12	GY-99540-57	
D. Amber narrow-mouth glass jugs				
2500 (80)	38-430	6	GY-99540-34	
4000 (135)	38-430	4	GY-99540-35	
E. Clear wide-mouth glass septa jars				
60 (2)	53 mm	24	GY-99540-20	
125 (4)	58 mm	24	GY-99540-21	

Cole-Parmer® Lab Markers

Fine-point fiber tip writes on cold, wet surfaces. Quick-drying ink permanently marks plastic; washes off glass, metal, or porcelain—won't bleed when subjected to water or alcohol.

Color	Qty/pk	Catalog number	Price
Black	10	GY-09964-00	
Blue		GY-09964-10	
Green		GY-09964-15	
Red		GY-09964-20	

**Precleaned Glass Sample Bottles****Save time—no need to prewash bottles**

- Vacuum and ionization process removes contaminants such as dirt, dust, carton lint, oil, aerosols, and fine glass particles
- Individually sealed for long-term storage with no contamination
- Bottles come with PTFE-lined caps; jars are available with PTFE-lined or pulp/vinyl-lined caps
- Clear wide-mouth square glass bottles have French square design for a more efficient use of shelf space while making the container easier to grip when pouring



Capacity mL (oz)	Neck size (mm)	Qty/ cs	Catalog number	Price/ cs
A. Clear narrow-mouth bottles				
30 (1)	20	48	GY-34606-00	
60 (2)	20		GY-34606-10	
120 (4)	22	24	GY-34606-20	
240 (8)	24		GY-34606-30	
480 (16)	28		GY-34606-40	
960 (32)	33	12	GY-34606-50	
B. Clear wide-mouth square glass bottles				
15 (½)	20		GY-34610-00	
30 (1)	24	48	GY-34610-10	
60 (2)	28		GY-34610-20	
120 (4)	33		GY-34610-30	
240 (8)	43	24	GY-34610-40	
480 (16)	48		GY-34610-50	
960 (32)†	58	12	GY-34610-60	
C. Amber narrow-mouth glass bottles				
60 (2)	20		GY-34607-00	
120 (4)	22	24	GY-34607-10	
240 (8)	24		GY-34607-20	
480 (16)	28		GY-34607-30	
960 (32)	33	12	GY-34607-40	
D. Amber wide-mouth glass bottles				
30 (1)	28		GY-34609-00	
60 (2)	33		GY-34609-10	
120 (4)	38	24	GY-34609-20	
240 (8)	45		GY-34609-30	
480 (16)	53	12	GY-34609-40	
960 (32)	53		GY-34609-50	
E. Clear wide-mouth glass jars with PTFE lined caps				
30 (1)	43	48	GY-34608-00	
60 (2)	53		GY-34608-10	
120 (4)	58	24	GY-34608-20	
270 (9)	70		GY-34608-30	
480 (16)	89	12	GY-34608-40	
960 (32)	89		GY-34608-50	
F. Clear wide-mouth glass jars with pulp/vinyl lined caps (not shown)				
30 (1)	43	48	GY-34608-05	
60 (2)	53		GY-34608-15	
120 (4)	58	24	GY-34608-25	
270 (9)	70		GY-34608-35	
480 (16)	89	12	GY-34608-45	
960 (32)	89		GY-34608-55	

†Graduated in mL and oz.

Wide-Mouth Sample Containers, LDPE

Store liquids or solids

- Include leakproof polypropylene screw caps



Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
10 (1/3)	22 mm	6	GY-62260-22	
30 (1)	22 mm		GY-62260-32	
90 (3)	35 mm		GY-62260-52	
125 (4)	54 mm		GY-62260-60	
180 (6)	54 mm		GY-62260-72	

Solution Bottles

Don't worry about stoppers fitting

- Tooled top for accurate stopper fit (stoppers not included)
- Meet ASTM E 438-92, Type 1, Class A requirements



Capacity L (gal.)	Stopper number	Qty/pk	Catalog number	Price/pk
1 (1/4)	6	6	GY-34500-00	
2 (1/2)	6	4	GY-34500-02	
5 (1 1/4)	10	1	GY-34500-04	
10 (2 1/2)	11.5	1	GY-34500-06	
20 (5 1/4)			GY-34500-08	

Heavy-Duty Solution Bottles

Eliminate cleaning up shattered glass

- KimCote™ plastic coating keeps broken glass in place for easy clean-up
- 5-gallon size meets Federal Spec DD-B-597
- Stoppers not included
- Autoclavable to 120°C (250°F).



Capacity L (gal.)	Stopper no.	Qty/pk	Catalog number	Price/pk
9.5 (2 1/2)	12	1	GY-34523-29	
13.2 (3 1/2)			GY-34523-31	
19 (5)			GY-34523-33	
45.5 (12)			GY-34523-35	

Safety-Coated Bottles

Non-slip grip for better handling

- Eliminate shattered glass with non-slip plastisol coating
- Resists acid and solvent degradation
- Meets USP Type III Standards
- Not autoclavable



Capacity L (oz)	Color	Cap size	Qty/cs	Catalog number	Price/cs
1/2 (16)	Clear	28–400	24	GY-08922-50	
1 (32)		33–400	12	GY-08922-60	
1/2 (16)	Amber	28–400	24	GY-08922-55	
1 (32)		33–400	12	GY-08922-65	
4 (128)		38–430	4	GY-08922-75	

Nalgene™ Square Wide-Mouth Bottles, HDPE

Maximize your sample storage space

- Space-saving design eliminates wasted shelf space
- Ideal for mailing
- Leakproof, PP screw enclosure



Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
60 (2)	28 mm	12	GY-06022-01	
175 (6)	38 mm		GY-06022-02	
250 (8)	43 mm		GY-06022-03	
500 (16)	53 mm		GY-06022-04	
1000 (32)	63 mm	6	GY-06022-05	
Replacement caps, PP				
—	28 mm	12	GY-63000-06	
—	43 mm		GY-63000-10	
—	53 mm		GY-63000-12	
—	63 mm		GY-63000-14	

Don't Forget...

Bottletop Dispensers
See pages 243–246



Bottles

Graduated Square Bottles, PC**Clearly see your samples**

- Know your remaining sample size with molded-in graduations
- Drip-proof pouring with spout design
- Autoclavable



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
125 (4) 250 (8)	38-430	6	GY-06040-50 GY-06040-55	
500 (16) 1000 (32)	38-430	4	GY-06040-60 GY-06040-65	
2000 (64)	53B	1	GY-06040-70	
Replacement caps, PP				
—	38-430	12	GY-63000-18	
—	53B		GY-63000-20	

Cole-Parmer® Graduated Rectangular Wide-Mouth Bottles, HDPE**Open with ease**

- Keep track of samples by attaching tag to closure ring
- Leakproof seal with included plastic insert
- Good chemical resistance



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
25 (¾) 50 (1¾) 100 (3½) 250 (8)	24 29 40 40	10	GY-06019-70 GY-06019-72 GY-06019-74 GY-06019-76	
500 (16) 1000 (32) 2000 (64)	53 66 66	6	GY-06019-78 GY-06019-80 GY-06019-82	

Square Wide-Mouth Bottles, PP**Fill and remove samples with ease**

- Leakproof, PP screw enclosures
- Withstand temperatures to 275°F (135°C)
- Autoclavable



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
60 (2) 175 (6) 250 (8)	28 mm 38 mm 43 mm	12	GY-06022-61 GY-06022-62 GY-06022-63	
500 (16) 1000 (32)	53 mm 63 mm	12 6	GY-06022-64 GY-06022-65	

Graduated Square Bottles, HDPE**Easy, drip-proof pouring**

- Know sample size with molded-in graduations
- Leakproof PP screw enclosures



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
30 (1) 60 (3) 125 (4)	20 mm 24 mm 38-430	12	GY-06252-00 GY-06252-10 GY-06252-20	
250 (8) 500 (16) 1000 (32)	38-430	12 6	GY-06252-30 GY-06252-40 GY-06252-50	
Replacement caps, PP				
—	24 mm	12	GY-63000-04	
—	38-430		GY-63000-18	

Wash Bottles, LDPE**No leaking while dispensing**

- One-piece polyallomer (PA) cap/delivery tube eliminates leaks
- Precise delivery with needle-thin stream



Capacity mL (oz)	Qty/ pk	Catalog number	Price/ pk
125 (4) 250 (8) 500 (16) 1000 (32)	12	GY-06053-21 GY-06053-41 GY-06053-61 GY-06053-81	

Multipurpose Jars, PC and PEI

Use as water baths or freeze-drying chambers

- Convert 06761-12, -25, and -30 to vacuum jars with optional vacuum plate
- Do not autoclave jars if using them for vacuum service



Capacity L (gal.)	Qty/ pk	Catalog number	Price/ pk
Polycarbonate jars, include lids			
2 (½) 4.75 (1¼) 8.5 (2¼) 19 (5)	1	GY-06761-20 GY-06761-12† GY-06761-00 GY-06761-27	
Amber polyetherimide jars, no lids			
8.25 (2½) 19 (5)	1	GY-06761-25† GY-06761-30†	

†These models fit the vacuum plates sold below.

[GY-06763-02](#) Vacuum plate for 06761-25, -30

[GY-06413-20](#) Tygon® vacuum tubing 1/4" ID. Pack of 10 ft (3.0 m)

Unitary™ Wash Bottles, LDPE

Fill up faster with wide-mouth design

- No need to shake or tip to dispense sample
- Increase flow by trimming tip
- One-piece leakproof design with PP screw closure



Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
125 (4)	24 mm	6	GY-06050-20	
250 (8)	38 mm	4	GY-06050-30	
500 (16)	38 mm	4	GY-06050-50	
750 (24)	38 mm	4	GY-06050-70	
1000 (32)	43 mm	2	GY-06050-90	

[GY-63000-04](#) Replacement caps, PP, 24 mm. Pack of 12

[GY-63000-10](#) Replacement caps, PP, 43 mm. Pack of 12

Wash Bottles, FEP

Won't contaminate reagents

- Resistant to most chemicals
- One-piece closure with spout prevents leaks
- Autoclavable



Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
125 (4)	24 mm	1	GY-06052-20	
250 (8)	24 mm		GY-06052-40	
500 (16)	38–430		GY-06052-60	
1000 (32)	38–430		GY-06052-80	

Economy Wash Bottles, LDPE

Eliminate pesky leaks

- One-piece molded PP closure and stem eliminate leaks
- Easily cut stem to increase flow



Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
125 (4)	24 mm	6	GY-62304-20	
250 (8)	24 mm		GY-62304-30	
500 (16)	28 mm		GY-62304-40	
1000 (32)	38–430	4	GY-62304-50	

“Right-to-Know” Safety Wash Bottles, LDPE

Labels will not wear off

- Nonmetallic vented closure prevents leaks
- Compliant with OSHA regulations



Labeled chemical	Closure color	Qty/ pk	Catalog number	Price/ pk
500-mL bottles, cap size 28 mm				
Acetone	Red	6	GY-62300-06	
Ethanol	White		GY-62300-16	
Methanol	Green		GY-62300-26	
Isopropanol	Yellow		GY-62300-36	
Distilled water	Natural		GY-62300-46	
500-mL assortment pack		5	GY-62300-56	
1000-mL bottles, cap size 38 mm				
Ethanol	White	4	GY-62300-17	
Methanol	Green		GY-62300-27	
Isopropanol	Yellow		GY-62300-37	
Distilled water	Natural		GY-62300-47	

Solvent Venting Wash Bottles, LDPE

Eliminate dripping and sputtering

- Fill bottles easily with wide-mouth design
- Easy identification with color-coding



Labeled chemical	Closure color	Qty/pk	Catalog number	Price/pk
250-mL bottles				
Acetone	Red	5	GY-62303-00	
Ethanol	Orange	5	GY-62303-10	
Isopropanol	Blue	5	GY-62303-20	
Methanol	Green	5	GY-62303-30	
Custom	Natural	5	GY-62303-80	
250-mL assortment pack		5	GY-62303-35	
500-mL bottles				
Acetone	Red	5	GY-62303-40	
Ethanol	Orange	5	GY-62303-50	
Isopropanol	Blue	5	GY-62303-60	
Methanol	Green	5	GY-62303-70	
Custom	Natural	5	GY-62303-90	
500-mL assortment (1 of each)		5	GY-62303-75	

Safety-Labeled Wash Bottles, LDPE

Keep solutions close at hand with easy identification

- Quickly identify contents by large lettering and color-coded caps
- These 500-mL bottles include a leakproof polyallomer (PA) cap and nozzle
- Assortment pack allows for easy ordering—includes one each of water, acetone, ethanol, isopropanol, methanol, and toluene†



Labeled chemical	Cap color	Qty/pk	Catalog number	Price/pk
Water	Blue	6	GY-06256-30	
Distilled water	Blue		GY-06256-20	
Deionized water	Blue		GY-06256-14	
Acetone	Red		GY-06256-40	
Methanol	Green	6	GY-06256-50	
Isopropanol	Yellow		GY-06256-60	
Toluene†	Red		GY-06256-70	
70% ethanol	Natural		GY-06256-13	
100% ethanol	Natural	6	GY-06256-81	
Bleach	Blue		GY-06256-90	
Methyl ethyl ketone	Green		GY-06256-10	
Soap	Blue		GY-06256-11	
Assortment pack for water, acetone, ethanol, isopropanol, methanol, and toluene†	—	6	GY-06256-00	

†Not recommended for long-term storage of toluene.

Cole-Parmer® Trigger Spray Bottle, HDPE

Get the stream you need

- Adjusts from fine mist to powerful stream
- Polypropylene cap and spray assembly with stainless steel bearings



Cole-Parmer®

Capacity mL (oz)	Qty/pk	Catalog number	Price/pk
650 (22)	4	GY-06091-05	

Dispenser Bottles, HDPE

Deliver a wide fine mist

- Easy-to-clean valve and intake
- Deliver thin liquids with ease



Capacity mL (oz)	Qty/pk	Catalog number	Price/pk
500 (16)	4	GY-06091-01	
1000 (32)		GY-06091-03	

Graduated Sample Containers, PP

No more leaks—leakproof for better sampling

- Easier sampling—graduated for measuring directly
- Durable design—translucent, unbreakable, and chemical resistant
- Sterile containers are individually wrapped



Capacity mL (oz)	Description	Qty/pk	Catalog number	Price/pk
118 (4)	Sterile, no tab	100	GY-06049-40	
118 (4)	Nonsterile with white tab	500	GY-06049-50	

Wide-Mouth Sample Containers, PS

Won't break like glass



- Get glass-like clarity with more rigidity
- Included screw caps keep samples secure
- Containers are polystyrene (PS); lids are polypropylene (PP) with polyethylene foam lines

Capacity mL (oz)	Dimensions (H x OD)		Cap size	Qty/ pk	Catalog number	Price/ pk
	in.	mm				
60 (2)	2 x 1¾	53 x 46	53 mm	48	GY-62263-20	
120 (4)	2¾ x 2	68 x 53	58 mm	36	GY-62263-30	
240 (8)	3½ x 2½	89 x 65	89 mm	36	GY-62263-40	
480 (16)	3¾ x 3½	96 x 89	89 mm	24	GY-62263-50	
960 (32)	4¾ x 3¾	120 x 98	120 mm	24	GY-62263-60	

Wide-Mouth Bottles, LDPE

Get better sample visibility



- Translucent material provides flexibility and impact resistance
- Ideal for shipping liquids—include leakproof polypropylene (PP) screw caps

Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
30 (1)	28 mm	12	GY-06043-10	
60 (2)	28 mm		GY-06043-20	
125 (4)	38 mm		GY-06043-30	
250 (8)	43 mm	12	GY-06043-40	
500 (16)	53 mm		GY-06043-50	
1000 (32)	63 mm	6	GY-06043-60	

Straight-Side, Wide-Mouth Containers, PC and PP

Stack easily to maximize space

- Easy viewing of samples for simple identification
- PP screw enclosures keep a tight seal
- Autoclavable (autoclave pieces separately)



Capacity mL (oz)	Cap size (mm)	Qty/ pk	Catalog number	Price/ pk
Polycarbonate (PC) containers				
30 (1)	43	4	GY-06760-50	
60 (2)	53		GY-06760-52	
125 (4)	70		GY-06760-54	
250 (8)	70	4	GY-06760-56	
500 (16)	120		GY-06760-58	
1000 (32)	120		GY-06760-60	
Polypropylene (PP) containers				
15 (½)	38	12	GY-06758-48	
30 (1)	43		GY-06758-50	
60 (2)	53		GY-06758-52	
125 (4)	70		GY-06758-54	
250 (8)	70	6	GY-06758-56	
500 (16)	120		GY-06758-58	
1000 (32)	120		GY-06758-60	

Cole-Parmer® Economical Bottle Carrier

Easily and safely transport solutions

- Lightweight, durable carrier transports up to six bottles measuring 3½" diameter or less



Cole-Parmer®

Hole diameter	Qty/ pk	Catalog number	Price/ pk
3½" (89 mm)	1	GY-06722-04	

Cole-Parmer® Lab Markers

Fine-point fiber tip writes on cold, wet surfaces. Quick-drying ink permanently marks plastic; washes off glass, metal, or porcelain—won't bleed when subjected to water or alcohol.

Color	Qty/pk	Catalog number	Price
Black	10	GY-09964-00	
Blue		GY-09964-10	
Green		GY-09964-15	
Red		GY-09964-20	



Bottles

Cole-Parmer® Wide-Mouth Bottles, LDPE, HDPE, PP

Wide-mouth opening provides easier filling and sampling

- Offered in pack quantities or bulk pack quantities to meet your specific requirements
- Reliable and sturdy
- Choose from three materials; all include a leakproof PP closure

HDPE Bottles are manufactured with FDA-compliant resin for noncontaminating storage of chemicals and food-grade products.

LDPE Bottles are translucent and ideal for food and general-use laboratory applications.

PP Bottles are rigid and offer superior chemical resistance; perfect for storage and sampling of liquids and solids. **Autoclavable.**



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk	Qty/ cs	Catalog number	Price/ cs
HDPE							
60 (2)	22 mm	12	GY-06045-18		72	GY-06045-23	
125 (4)	28 mm	12	GY-06045-24		72	GY-06045-25	
250 (8)	28 mm	12	GY-06045-26		72	GY-06045-27	
500 (16)	28 mm	12	GY-06045-28		48	GY-06045-29	
1000 (32)	28 mm	6	GY-06045-44		24	GY-06045-31	
2000 (64)	45 mm	1	GY-06045-46		6	GY-06045-33	
LDPE							
60 (2)	22 mm	12	GY-06045-02		72	GY-06045-03	
125 (4)	28 mm	12	GY-06045-04		72	GY-06045-05	
250 (8)	28 mm	12	GY-06045-06		72	GY-06045-07	
500 (16)	28 mm	12	GY-06045-08		48	GY-06045-09	
1000 (32)	28 mm	6	GY-06045-14		6	GY-06045-11	
PP							
30 (1)	22 mm	12	GY-06045-48		72	GY-06045-41	
60 (2)	22 mm	12	GY-06045-56		72	GY-06045-43	
125 (4)	28 mm	12	GY-06045-58		72	GY-06045-45	
250 (8)	28 mm	12	GY-06045-62		72	GY-06045-47	
500 (16)	28 mm	12	GY-06045-66		48	GY-06045-49	
1000 (32)	28 mm	6	GY-06045-68		24	GY-06045-51	
2000 (64)	45 mm	1	GY-06045-74		6	GY-06045-53	

Sample Bottles, HDPE

Store samples without breaking your budget

- Provide good chemical resistance for a wide range of uses
- Bottles are made of HDPE; included cap is made of polypropylene and lined with polyethylene foam
- Easy mailing with optional mailing tube (sold separately)



Capacity mL (oz)	Qty/ pk	Catalog number	Price/ pk
180 (6)	300	GY-06043-01	
250 (8)	329	GY-06049-01	
500 (16)	198	GY-06050-05	

[GY-06049-30](#) Mailing tubes for bottle 06049-01. Carton of 240

Amber Wide-Mouth Bottles, HDPE

Keep light-sensitive samples safe

- Reduce UV light transmission
- Meet USP requirements
- Leakproof, PP closures



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
30 (1)	28 mm	12	GY-06027-50	
60 (2)	28 mm		GY-06027-52	
125 (4)	38 mm		GY-06027-54	
250 (8)	43 mm		GY-06027-56	
500 (16)	53 mm		GY-06027-58	
1000 (32)	63 mm	6	GY-06027-60	

Economy Wide-Mouth Bottles HDPE

Rugged design fits all your needs

- Leakproof, PP enclosure

Capacity mL (oz)	Cap size	Qty/pk	Cat. no.	Price/pk	Qty/cs	Cat. no.	Price/cs
30 (1)	28 mm	12	GY-06035-31		72	GY-06035-33	
60 (2)	28 mm		GY-06035-35			GY-06035-37	
125 (4)	38 mm		GY-06035-39			GY-06035-41	
250 (8)	43 mm	12	GY-06035-43		72	GY-06035-45	
500 (16)	53 mm	6	GY-06035-47		48	GY-06035-49	
1000 (32)	63 mm	6	GY-06035-51		24	GY-06035-53	



Cole-Parmer® Cylindrical Jars, HDPE

Get easy sample access

- Wide neck for easy powder sampling
- Include screw caps to keep samples intact
- Won't break—made of high-density polyethylene for sturdy storage

Capacity mL (oz)	Dimensions (H x OD)		Qty/ pk	Catalog number	Price/pk
	in.	mm			
70 (2.4)	2 ⁵ / ₁₆ x 2	60 x 50	10	GY-63100-10	
120 (4)	2 ³ / ₄ x 2 ¹ / ₄	71 x 56		GY-63100-12	
250 (8)	3 ³ / ₄ x 2 ³ / ₄	94 x 69		GY-63100-14	
500 (16)	4 ¹ / ₄ x 3 ³ / ₈	107 x 86		GY-63100-16	
1000 (32)	5 x 4 ³ / ₈	128 x 111	10	GY-63100-18	
1500 (51)	7 ³ / ₁₆ x 4 ³ / ₈	182 x 111		GY-63100-20	
2000 (64)	9 ¹ / ₄ x 4 ³ / ₈	235 x 111		GY-63100-22	



4-L (1-Gal.) Square Bottles, HDPE and PP

Get a better grip

- Molded-in handgrips for safer handling
- Suitable for solids and other specimens
- HDPE is leakproof; PP is autoclavable



Wide-Mouth Bottles, PP

Ideal for sensitive sample requirements

- Autoclavable bottles provide exceptional sample protection
- Include leakproof, polypropylene (PP) screw caps



Material	Capacity L (gal.)	Cap size	Qty/ pk	Catalog number	Price/ pk
HDPE	4 (1)	100 mm	1	GY-06022-97	
PP				GY-06022-95	
Replacement caps, PP					
—	—	100 mm	12	GY-63000-16	

Capacity, mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
30 (1)	28 mm	12	GY-06048-01	
60 (2)	28 mm		GY-06048-02	
125 (4)	38 mm		GY-06048-10	
250 (8)	43 mm	6	GY-06048-20	
500 (16)	53 mm		GY-06048-30	
1000 (32)	63 mm		GY-06048-40	

Accessories

- [GY-63000-06](#) Replacement cap, PP, 28 mm. Case of 12
- [GY-63000-10](#) Replacement cap, PP, 43 mm. Case of 12
- [GY-63000-12](#) Replacement cap, PP, 53 mm. Case of 12
- [GY-63000-14](#) Replacement cap, PP, 63 mm. Case of 12

Economy Amber Wide-Mouth Bottles, HDPE

Easy sampling of light-sensitive materials

- Leakproof, PP closure
- Keep light-sensitive materials safe



Wait!

**There's More at
[ColeParmer.com](#)**

Find more than
100,000 products
to choose from!

Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk	Qty/ cs	Catalog number	Price/ cs
30 (1)	28 mm	12	GY-06036-01		72	GY-06036-03	
60 (2)	28 mm		GY-06036-05		72	GY-06036-07	
125 (4)	38 mm		GY-06036-09		72	GY-06036-11	
250 (8)	43 mm	12	GY-06036-13		72	GY-06036-15	
500 (16)	53 mm	6	GY-06036-17		48	GY-06036-19	
1000 (32)	63 mm	6	GY-06036-21		24	GY-06036-23	

Bottles

**Wide-Mouth
PassPort™ IP2
Bottles, HDPE**

Ship hazardous
samples reliably

- Safe for small quantities of corrosives or cleaning supplies
- Leakproof, PP screw cap
- Volume, registration notch, material, and IP2 markings molded into bottom



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
125 (4)	38 mm	12	GY-62150-00	
250 (8)	43 mm		GY-62150-10	
500 (16)	53 mm		GY-62150-20	
1000 (32)	63 mm	6	GY-62150-30	

Wide-Mouth Jars, HDPE

Get better access to
your sample

- Wide mouth design is perfect for powder sampling
- Include HDPE closure



Capacity L (gal.)	Neck dia	Qty/ pk	Catalog number	Price/ pk
5 (1 1/3)	110 mm	1	GY-06046-10	
10 (2 2/3)			GY-06046-20	

**Wide-Mouth Bottles,
HDPE and PP**

Sample easier with
wider mouth

- Store dry chemicals or take water samples quickly
- Include PP caps




Capacity L (gal.)	Cap size	Qty/ pk	Catalog number	Price/ pk
High-density polyethylene (HDPE) bottles, leakproof				
2 (½)	100 mm	1	GY-06055-10	
4 (1)	100 mm		GY-06055-20	
Polypropylene (PP) bottles. Autoclavable but not leakproof				
2 (½)	100 mm	1	GY-06057-10	
4 (1)	100 mm		GY-06057-20	
Replacement PP cap. Autoclavable				
—	100 mm	12	GY-63000-16	

**Fluorinated Wide-Mouth
Bottles, FLPE**

Stop worrying about
splitting and puncturing

- Improved barrier for reduced solvent absorption and penetration
- Leakproof FLPP closures
- Stand up to many solvents



Fluorinated
bottles are
ideal for
many
solvents.



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
125 (4)	38 mm	12	GY-62151-00	
250 (8)	43 mm		GY-62151-10	
500 (16)	53 mm		GY-62151-20	
1000 (32)	63 mm	6	GY-62151-30	
2000 (64)	100 mm	1	GY-62503-00	

Wide-Mouth Jars, PET

Maintain
FDA-compliance

- Unbreakable design
- Resist moisture and oxygen to keep samples fresh



“The little jar that works!”



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
180 (6)	53 mm	12	GY-06043-70	
240 (8)	53 mm		GY-06043-71	
350 (12)	63 mm		GY-06043-72	
480 (16)	63 mm		GY-06043-73	
960 (32)	89 mm	12	GY-06043-74	
2000 (64)	100 mm	6	GY-06043-75	
4000 (128)	100 mm	4	GY-06043-76	

Wide-Mouth Bottles, FEP

Maintain sterility of samples

- Inert, FEP bottles withstand repeated sterilization
- Use in broad temperature range –157 to 302°F (–105 to 150°C)
- Autoclavable



Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
125 (4)	33 mm	1	GY-06044-10	
250 (8)	43 mm		GY-06044-20	
500 (16)	48 mm		GY-06044-30	
1000 (32)	53 mm		GY-06044-40	

Economy Wide-Mouth Bottles, PPCO

Save money while keeping samples contained

- Lightweight design minimizes added weight
- Leakproof PP closures
- Autoclavable.




Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk	Qty/cs	Catalog number	Price/cs
30 (1)	28 mm	12	GY-06034-51		72	GY-06034-53	
60 (2)	28 mm		GY-06034-55			GY-06034-57	
125 (4)	38 mm		GY-06034-59			GY-06034-61	
250 (8)	43 mm		GY-06034-63			GY-06034-65	
500 (16)	53 mm	6	GY-06034-67		48	GY-06034-69	
1000 (32)	63 mm	6	GY-06034-71		24	GY-06034-73	

Cole-Parmer® Graduated Bottles, PP

Don't ever guess what's in a bottle again

- Large write-on area provides easy labeling
- Know how much sample remains with printed graduations
- Chemical resistant and leakproof
- Autoclavable



Cole-Parmer® Wide-Mouth Bottles, PTFE

Store corrosive materials without worry

- Withstand temperatures up to 540°F (280°C)
- Leakproof design
- Autoclavable



Cole-Parmer®

Cole-Parmer®

Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
60 (2)	33 mm	12	GY-06041-10	
125 (4)	38 mm		GY-06041-12	
250 (8)	45 mm		GY-06041-14	
500 (16)	55 mm		GY-06041-16	
1000 (32)	65 mm	6	GY-06041-18	

Capacity mL (oz)	Cap size	Qty/pk	Catalog number	Price/pk
1 (0.03)	10	1	GY-06300-12	
5 (1.6)	10		GY-06300-16	
10 (0.3)	12		GY-06300-18	
25 (0.8)	19	1	GY-06300-24	
50 (1.6)	25		GY-06300-26	
100 (3)	35		GY-06300-28	
150 (5)	35		GY-06300-34	
250 (8)	42	1	GY-06300-36	
500 (16)	52		GY-06300-38	
1000 (32)	57		GY-06300-42	
2000 (64)	60		GY-06300-44	

Cole-Parmer® Lab Markers

Fine-point fiber tip writes on cold, wet surfaces. Quick-drying ink permanently marks plastic; washes off glass, metal, or porcelain—won't bleed when subjected to water or alcohol.

Color	Qty/pk	Catalog number	Price
Black	10	GY-09964-00	
Blue		GY-09964-10	
Green		GY-09964-15	
Red		GY-09964-20	



Bottles

Mason Containers, PP

Store samples safely
for extended periods

- Heavy-walled for added strength
- Included PP screw caps



Capacity mL (oz)	Cap size	Qty/ pk	Catalog number	Price/ pk
480 (16)	70 mm	6	GY-06029-11	
960 (32)		6	GY-06029-21	
2000 (64)		3	GY-06029-31	
4000 (128)		2	GY-06029-41	
8000 (256)		1	GY-06029-50	

Glass Sample Jars

No precleaning required

- Precleaned for semivolatiles, pesticides, PCBs, and metals
- Include chemical-resistant PTFE-lined screw cap



Capacity mL (oz)	Neck size	Qty/ cs	Catalog number	Price/ cs
120 (4)	53 mm	24	GY-34605-30	
240 (8)	70 mm		GY-34605-55	
480 (16)	89 mm	12	GY-34605-60	
960 (32)	89 mm		GY-34605-70	

Wide-Mouth Bottles

Maintain sample integrity

- Tamper-evident labels
- Ideal for soil samples and specimens
- Included caps



Capacity mL (oz)	Cap material	Clear bottles				Amber bottles			
		Cap size	Qty/cs	Cat. no.	Price/cs	Cap size	Qty/cs	Cat. no.	Price/cs
30 (1)	Foamed poly liner	33–400	48	GY-34630-51		28–400	24	GY-34632-31	
60 (2)		38–400	48	GY-34630-52		33–400	24	GY-34632-32	
125 (4)		48–400	24	GY-34630-53		38–400	24	GY-34632-33	
250 (8)		58–400	24	GY-34630-54		45–400	24	GY-34632-34	
500 (16)		70–400	24	GY-34630-55		53–400	12	GY-34632-35	
1000 (32)		—	—	—		53–400	12	GY-34632-36	
30 (1)	Polyvinyl liner	33–400	48	GY-34630-61		28–400	24	GY-34632-41	
60 (2)		38–400	48	GY-34630-62		33–400	24	GY-34632-42	
125 (4)		48–400	24	GY-34630-63		38–400	24	GY-34632-43	
250 (8)		58–400	24	GY-34630-64		45–400	24	GY-34632-44	
500 (16)		70–400	24	GY-34630-65		53–400	12	GY-34632-45	
1000 (32)		—	—	—		53–400	12	GY-34632-46	

Graduated Bottles

Less spilling during
transfer

- Wide-mouth design with PTFE-lined caps
- Easily measure samples as you sample and pour—graduations in both mL and ounces



Capacity, mL (oz)	Graduations, mL (oz)	Neck size	Qty/cs	Catalog number	Price/cs
30 (1)	5 (1/4)	33 mm	48	GY-34601-10	
60 (2)		38 mm	48	GY-34601-20	
120 (4)		48 mm	24	GY-34601-30	
240 (8)	20 (1/4)	58 mm	24	GY-34601-40	
480 (16)		70 mm	24	GY-34601-50	
960 (32)		58 mm	12	GY-34601-60	

Safety-Coated Sample Jars

Will not shatter

- Reduce breakage with safety-coated design
- Chemical-resistant, PTFE-lined caps for tight seal



Style	Capacity mL (oz)	Neck size	Qty/ cs	Catalog number	Price/ cs
Straight sided	260 (8)	70 mm	24	GY-34604-10	
	480 (16)		24	GY-34604-20	
	960 (32)		12	GY-34604-30	
Standard rounded	1920 (64)	83 mm	6	GY-34604-40	