

Beaker and Bottle Brushes

No more dirty glassware—clean it easily

- Choose the right brush for your glassware from shape to bristle type
- Radial tip brushes reduce scratches and provide effective cleaning
- Contoured brushes allow you to reach angled areas with ease
- Tufted end brushes fold to get in narrow-mouth vessels and then open once inside
- Choose angle and swivel head brushes to easily clean curves
- Volumetric models have a bendable handle to fit container's shape

Wait!

**There's More at
ColeParmer.com**

Find more than
100,000 products
to choose from!



Key	Description	Brush bristles			Handle length		Catalog number	Qty/pk	Price
		Size (L x W)		Material	in.	cm			
Beaker/bottle/flask brushes									
A	Radial tip, small	3" x 1¼"	7.6 x 3.2 cm	Polyester	9	22.9	GY-84550-80	12	
	Radial tip, large	6" x 3⅞"	15.2 x 7.9 cm		10	25.4	GY-84550-81	12	
	Radial tip, XL	6" x 3⅞"	15.2 x 7.9 cm		13	33	GY-84550-82	6	
B	Contoured, large	10" x 3⅞"	25.4 x 7.9 cm	Natural	7	17.8	GY-84550-83	6	
C	Fan tip, small	6" x 1½"	15.2 x 3.8 cm		6	15.2	GY-84550-86	12	
D	Tufted end, small	4¼" x 1½"	11.4 x 3.8 cm		7½	19.1	GY-84550-87	12	
	Tufted end, XL	4¼" x 1¾"	11.4 x 4.5 cm	19½	49.5	GY-84550-88	12		
E	Angled head, small	3" x 1"	7.6 x 2.5 cm	Polyester	13	33	GY-84550-76	1	
	Angled head, large	4" x 1½"	0.2 x 3.8 cm	Natural	12	30.5	GY-84550-89	12	
	Volumetric, large	4" x 1½"	10.2 x 3.8 cm	Nylon	10	25.4	GY-84550-90	12	
Bottle/cylinder brushes									
F	Radial tip, white bristles	6½" x 3⅞" dia	16.5 x 7.9 cm dia	Nylon	9	22.9	GY-84550-00	1	
	Radial tip, X-shaped yellow bristles						GY-84550-03	1	
	Radial tip, X-shaped green bristles						GY-84550-04	1	
G	Radial tip, white bristles	6½" x 3⅞" dia	16.5 x 7.9 cm dia	Polyester	9	22.9	GY-84550-70	1	
	Radial tip, red bristles						GY-84550-72	1	
	Radial tip, blue bristles						GY-84550-74	1	
H	Radial tip, small	4½" x 1¼"	11.4 x 4.5 cm	Natural	7	17.8	GY-84550-15	12	
	Radial tip, medium	5" x 2"	12.7 x 5.1 cm	Nylon	11	27.9	GY-84550-16	12	
	Radial tip, large	6" x 2⅝"	15.2 x 6.1 cm		14	35.6	GY-84550-17	12	
	Radial tip, XL	6" x 2¾"	15.2 x 7 cm		18	45.7	GY-84550-18	12	
	Radial tip, XXL				24	60.9	GY-84550-19	12	
Glassware brushes									
I	Foam-tipped	7" x 2"	17.8 x 5.1 cm	Polypropylene/foam	13	33	GY-84552-14	1	
		5" x 2"	12.7 x 5.1 cm		7	17.8	GY-84552-16	1	
J	Cylinder/Funnel brush	16" x 3⅞"	40.6 x 7.9 cm	Polyester	10½	26.7	GY-84550-71	12	
		24" x 3⅞"	60.9 x 7.9 cm		17½	44.5	GY-84550-73	12	
K	Swivel head brush	1¾" dia	3.2 cm	Natural	13	33	GY-84550-55	1	

Test Tube and Pipette Brushes

Get even the smallest diameter vessel clean

- Long, flexible handles to fit in narrow areas
- Clean all sides with ease by using radial tip brushes with sprays on tip and sides
- Small openings are easy to clean with tied tip brushes—expand to clean sides and bottom without scratching
- Get a more thorough cleaning with sponge tips—expand for better coverage
- Fit all tube sizes with adjustable test tube brush—spreads apart to clean all sides



Key	Description	Brush bristles		Material	Handle length		Catalog number	Qty/pk	Price
		Size (L x W)			in.	cm			
Test tube brushes									
A	Radial tip	4" x ¾"	10.2 x 1.9 cm	Nylon	6	15.2	GY-84550-52	12	
		3" x ½"	7.6 x 1.3 cm		5	12.7	GY-84551-50	12	
		3" x 1"	7.6 x 2.5 cm		7	17.8	GY-84551-51	12	
		4" x 1⅝"	10.2 x 4.1 cm	Natural	9	22.9	GY-84551-52	12	
		3" x ⅝"	7.6 x 1.6 cm		9	22.9	GY-84551-53	12	
		3" x ⅞"	7.6 x 2.2 cm		9	22.9	GY-84551-54	12	
	Radial tip/SS wire stem	4" x 1⅝"	10.2 x 3.5 cm	Nylon	7	17.8	GY-84551-55	12	
		3" x ½"	7.6 x 1.3 cm		5	12.7	GY-84551-56	12	
		3½" x ¾"	8.9 x 1.9 cm		5½	14	GY-84551-57	12	
		4" x 1⅝"	10.2 x 3.5 cm		7	17.8	GY-84551-58	12	
Sponge tip	3" x 1⅝"	7.6 x 3.5 cm	Natural	6	15.2	GY-84551-59	12		
	3" x 1½"	7.6 x 3.8 cm		7	17.8	GY-84551-61	12		
B	Tied tip	3½" x ½"	8.9 x 1.3 cm	Natural	7	17.8	GY-84551-62	12	
		6" x ¾"	15.2 x 1.9 cm		6	15.2	GY-84551-63	12	
		3" x 1¼"	7.6 x 3.2 cm		9	22.9	GY-84551-64	12	
C	Adjustable	6¾" x ½"	17.1 x 1.3 cm	Nylon	2¼	5.7	GY-84551-65	1	
Centrifuge/micro brushes									
D	Tied tip	3" x ⅜ to ½"	7.6 x 0.5 to 1.3 cm	Natural/nylon	6	15.2	GY-84551-90	12	
E	Radial tip	3½" x ⅜ to 1½"	8.9 x 0.5 to 1.3 cm	Natural	5½	14	GY-84551-91	12	
F		3" x ½ to 1⅜"	7.6 x 1.0 to 0.5 cm		6	15.2	GY-84551-92	12	
F	Tied tip	4½" x ¾ to 1½"	11.4 x 1.9 to 3.8 cm		7½	19.1	GY-84551-93	12	
G	Micro tube	3" x ¼"	7.6 x 0.6 cm	Nylon	3	7.6	GY-84551-94	12	
	Micro tube, fan tip	3" x ⅜ to ½"	7.6 x 0.5 to 1.3 cm		4	10.2	GY-84551-95	12	
Pipette/burette brushes									
H	Tapered tip	5" x ⅝ to ⅞"	12.7 x 1.6 to 2.2 cm	Natural	31	78.7	GY-84550-56	12	
I		6" x ⅜"	15.2 x 0.5 cm	Nylon	9	22.9	GY-84550-60	12	
		7½" x ⅙ to ¾"	19.1 x 0.3 to 1.9 cm		9	22.9	GY-84550-61	12	
		7½" x ⅜ to ¾"	19.1 x 0.5 to 1.9 cm		9½	24.1	GY-84550-62	12	
		7½" x ⅜ to ⅝"	19.1 x 0.5 to 0.8 cm		9½	24.1	GY-84550-63	12	
		1½" x ⅙ tip; 4" x ¾" body	3.8 x 0.3 cm tip; 10.2 x 1.9 cm body		12	30.5	GY-84550-64	12	

Wait!

There's More at
ColeParmer.com

Browse **100,000+** products plus find
Technical Resources to help you choose