

Carboys, PC

Maintain visibility of your solution

- Know your fluid levels using 1-gallon and 5-liter increments
- Transport easily with leakproof, polypropylene cap
- Easily dispense with spigot models
- Autoclavable



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
Carboys				
10 (2½)	83B	GY-06063-60	1	
20 (5)	83B	GY-06063-62		
Carboys with spigot				
10 (2½)	83B	GY-06063-80	1	
20 (5)	83B	GY-06063-82		

Accessories

[GY-06070-22](#) Replacement spigot, PP

[GY-06258-30](#) Replacement caps; PP, 83B. Pack of 2

Cole-Parmer® Economical Heavy-Duty Carboys with Shoulder Handle, PP

Sturdy, octagonal design for added rigidity

- Easier filling and pouring with offset neck
- Leakproof spigot models have tube fitting for easier filling of containers
- Polypropylene closure reduces evaporation
- Autoclavable



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
Heavy-duty carboys with shoulder handle				
5 (1¼)	45 mm	GY-62507-00	1	
10 (2½)	90 mm	GY-62507-10	1	
20 (5¼)	90 mm	GY-62507-20	1	
Heavy-duty carboys with shoulder handle and spigot				
5 (1¼)	45 mm	GY-62507-05	1	
10 (2½)	90 mm	GY-62507-15	1	
20 (5¼)	90 mm	GY-62507-25	1	

Round Carboys, LDPE

No more leaks

- Leakproof, polypropylene cap for no worry transport
- Make transport easy with easy-grip handles
- Withstand corrosive materials
- Maintain adequate levels of solution—larger than 10-L carboys indicate 1-gallon and 5-liter increments
- Dispense easily with spigot models



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
Carboys with shoulder handles				
10 (2½)	83B	GY-06031-32	1	
15 (4)	83B	GY-06031-42		
20 (5)	83B	GY-06031-52		
25 (6½)	83B	GY-06031-62		
50 (13)	83B	GY-06031-72		
Wide-mouth carboys with shoulder handles				
10 (2½)	100 mm	GY-06056-22	1	
15 (4)	100 mm	GY-06056-32		
20 (5)	100 mm	GY-06056-52		
Carboy with spigot and shoulder handles				
4 (1)	38-430	GY-06065-10	1	
10 (2½)	83B	GY-06065-22		
20 (5)	83B	GY-06065-32		
25 (6½)	83B	GY-06065-52		
50 (13)	83B	GY-06065-62		
Carboys with ½" OD tubulation and shoulder handles				
10 (2½)	83B	GY-06080-42	1	
20 (5¼)	83B	GY-06080-52		

Accessories

[GY-06258-30](#) Replacement cap, PP, 83B. Pack of 2

[GY-63000-16](#) Replacement cap, PP, 100 mm. Pack of 12

[GY-63000-18](#) Replacement cap, PP, 38-430. Pack of 12

[GY-06070-22](#) Replacement spigot, PP

Carboys

Carboys, PP

Don't compromise sterility—
completely autoclavable

- Eliminate contaminations—select C-Flex® or Sani-Tech® high-purity closures
- No more difficult transporting—handles on top and bottom make carrying easy
- No more leaks—EZ Top® closure provides leakproof top connections
- Chemical resistant polypropylene (PP) is BPA-free and can be heated to 275°F (135°C)

NEW

Liters (gallons)	Cap size	Spigot	Catalog number	Price/pk
10 (2.6)	83 mm	No	GY-62500-01	
		Yes	GY-62500-03	
20 (5.3)	120 mm	No	GY-62500-05	
		Yes	GY-62500-07	
40 (10.5)	120 mm	No	GY-62500-09	
		Yes	GY-62500-11	

High-Purity EZ Top® Container Closures

Fits carboy size	Material	No. of ports	Catalog number	Price
10 L (2.6 gal.)	C-Flex	2	GY-62500-47	
	Sani-Tech	3	GY-62500-53	
20 L (5.3 gal.)	C-Flex	2	GY-62500-49	
	Sani-Tech	3	GY-62500-55	
40 L (10.5 gal.)	C-Flex	2	GY-62500-51	
	Sani-Tech	3	GY-62500-57	

[GY-62500-41](#) Replacement EZ Top closed top cap, 83 mm. 2/pk
[GY-62500-43](#) Replacement EZ Top closed top cap, 120 mm. 2/pk
[GY-62500-45](#) Replacement spigot, fits all sizes. 2/pk

Dispensing Jug
with Handle, HDPE

Save space

- Maximize storage space with unique cube design
- Keep aqueous solutions close at hand with easy dispensing
- Know sample levels with molded-in graduations
- Threaded LDPE spigot features a 35-mm stem for easy dispensing

NEW

Liters (gallons)	Cap size	Catalog number	Price/pk
20 (5)	70 mm	GY-06034-00	

FDA-Compliant Carboys, PP

Get a better grip with
two-handle construction

- Dispense solutions more quickly with optional spigot models—accept 1/4", 1/2", or 5/8" ID tubing
- Leakproof polypropylene (PP) cap
- Autoclavable



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
FDA-compliant carboys, USP Class VI				
10 (2½)	83B	GY-06063-12	1	
20 (5)	83B	GY-06063-22		
50 (13)	83B	GY-06063-32		
FDA-compliant carboys with spigot, USP Class VI				
10 (2½)	83B	GY-06064-50	1	
20 (5)	83B	GY-06064-60		
50 (13)	83B	GY-06064-70		

Cole-Parmer® Heavy-Walled
Carboys, HDPEDon't worry about
shipping!

- Resistant to most chemicals
- Effortless transport via heavy-duty swing handle

Narrow-Mouth Carboys
 come with rubber O-ring and HDPE closure with tag holes for identification. Available with or without leakproof spigot.

Wide-Mouth Carboys
 feature LDPE neck insert plug for maximum security.



Cole-Parmer®

Capacity L (gal.)	Cap size mm	Catalog number	Qty/pk	Price/pk
Narrow-mouth carboy				
5 (1¼)	46	GY-06035-93	1	
10 (2½)	55	GY-06035-95	1	
25 (6½)	55	GY-06035-97	1	
50 (13)	55	GY-06035-99	1	
Narrow-mouth carboys with spigot				
5 (1¼)	46	GY-06035-81	1	
10 (2½)	55	GY-06035-83	1	
25 (6½)	55	GY-06035-85	1	
50 (13)	55	GY-06035-87	1	
Wide-mouth carboys with neck insert				
5 (1¼)	90	GY-06035-89	1	
10 (2½)	120	GY-06035-91	1	

[GY-06035-76](#) Replacement spigot

Cole-Parmer® Economical Square Carboys, HDPE

Handle and pour contents easily

- Molded-in handle and indented bottom for easy grabbing
- No leaks with leakproof HDPE threaded cap
- Get easier dispensing by drilling and attaching optional faucet
- Autoclavable; withstand temperature from 68 to 248°F (20 to 120°C)



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
10 (2½)	63 mm	GY-06031-80	1	
20 (5)		GY-06031-82		
10 (2½)	63 mm	GY-06031-84	4	
20 (5)		GY-06031-86		

Accessories

[GY-06031-88](#) HDPE faucet, ¾"

Vacuum Carboys with Handles, PP

Withstand vacuum for up to eight hours

- Save money compared to stainless steel and fiberglass-reinforced containers
- No leaks with polypropylene leakproof cap
- Autoclavable



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
10 (2½)	83B	GY-06257-50	1	
20 (5)		GY-06257-60		

Accessories

[GY-06258-30](#) Replacement cap; PP, 83B. Pack of 2

Jerry Cans, HDPE/Quoral®

More environmentally friendly design

- Great alternative to glass, metal, or fluorinated containers
- Safely transport chemicals, solvents, adhesives, and cleaning solutions
- Use in temperatures from –148 to 248° (–100 to 120°C)



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
5 (1⅓)	53 mm	GY-06069-01	1	
10 (2½)	53 mm	GY-06069-03	1	
20 (5¼)	63 mm	GY-06069-05	1	

Safety Dispensing Jugs, LDPE

Just tip and pour

- Get a better grip with top and bottom handles for easier pouring
- Leakproof design for better transport
- Don't misplace spout plug—attached for easier locating



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
10 (2½)	83B	GY-06077-22	1	
20 (5)		GY-06077-52		

Accessories

[GY-06258-30](#) Replacement cap; PP, 83B. Pack of 2

Carboys

Rectangular Carboys, PP**Maintain sample integrity—autoclavable**

- No dead space—rectangular design does not waste space
- Leakproof for better handling
- Easy transport via steel handle
- Withstand temperatures to 248°F (120°C)



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
Rectangular carboys				
9 (2)	100 mm	GY-06082-60	1	
20 (5)	100 mm	GY-06082-70	1	
Rectangular carboys with spigot				
9 (2)	100 mm	GY-06082-80	1	
20 (5)	100 mm	GY-06082-90	1	

[GY-06070-22](#) Replacement spigot, PP

[GY-63000-16](#) Replacement cap; PP, 100 mm. Pack of 12

Rectangular Carboys, PC**Don't lose sight of your solution**

- Transparent, polycarbonate (PC) carboys allow easy viewing of solutions
- Easily transport carboys with steel handle
- Dispense samples quickly with spigot models
- Autoclavable



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
Rectangular carboy				
20 (5)	100 mm	GY-06084-70	1	
Rectangular carboys with spigot				
9 (2)	100 mm	GY-06084-80	1	
20 (5)	100 mm	GY-06084-90	1	

[GY-06070-22](#) Replacement spigot, PP

[GY-63000-16](#) Replacement cap; PP, 100 mm. Pack of 12

Rectangular Carboys, HDPE**Don't waste valuable bench space**

- Space-saving design for no dead area like traditional round carboys
- Leakproof high-density polyethylene (HDPE) for safe transport
- Liter/gallon graduations for easy volume determination

Carboys with quick-action spigots accept $\frac{5}{8}$ " ID tubing, come with adapters for $\frac{1}{4}$ " to $\frac{1}{2}$ " ID tubing, and feature an integrally molded threaded boss. Aspirator carboys feature integral tubulation for $\frac{1}{2}$ " ID tubing.



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
Rectangular carboys				
9 (2)	100 mm	GY-06032-20	1	
20 (5)	100 mm	GY-06032-50	1	
Rectangular carboys with spigot				
9 (2)	100 mm	GY-06066-20	1	
20 (5)	100 mm	GY-06066-50	1	
Rectangular carboys with tubulation				
20 (5)	100 mm	GY-06082-50	1	

Filling/Venting Caps, PP**Quick and easy aseptic liquid transfer**

Filling/Venting Caps have thermoplastic elastomer (TPE) gaskets to allow aseptic liquid transfer. Include 6" (15 cm) and 30" (76 cm) lengths of platinum-cured silicone tubing. Autoclavable.

Fitting Kits include two PP fittings, two acetal nuts, two silicone gaskets, thermoplastic elastomer (TPE) port cap, instructions, and template. Autoclavable.



Tubing ID	Ports	Cap size	Catalog number	Qty/pk	Price/pk
Filling/venting caps					
½"	2 tubing, 1 venting	83B	GY-06258-00	1	
¼"		53B	GY-06258-10	1	
¼"		83B	GY-06258-20	1	
Barbed bulkhead fitting kits					
¼"	—	—	GY-06259-10	2	
½"	—	—	GY-06259-00	2	

Round Fluorinated Carboys with Grips, FLPE

Minimize solvent absorption and penetration

- Won't leak or puncture
- Know solution levels—graduated in 1-gallon and 5L increments
- Include FLPP cap



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
10 (2½)	83B	GY-62500-50	1	
20 (5)		GY-62500-60		

[GY-63000-20](#) Replacement cap; PP, 53B. Pack of 12

Lowboys™, HDPE and PP

Space-saving design is ideal for shelves

- Removable spigot provides quick dispensing
- Easily attach tubing to spigot with included HDPE tubing adapter



Tubing adapter (HDPE) included

Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
Lowboys, HDPE				
8 (2)	63 mm	GY-06059-05	1	
15 (4)		GY-06059-10		
Lowboys, PP, Autoclavable—remove spigot before autoclaving				
8 (2)	63 mm	GY-06059-25	1	
15 (4)		GY-06059-30		

[GY-06070-22](#) Replacement spigot, PP

Modular Carboy Racks

Keep fluid supplies close at hand

- Assemble racks in minutes without tools
- Heat and chemical resistant, easy to maintain, and rust free
- Made from a polypropylene blend for strength and rigidity

Each shelf can hold a five-gallon carboy as shown. Maximum temperature is 275°F (135°C).

Carboys 06060-10 shown on racks; order separately below.



Shelves	Overall dimensions (W x H x D)		Catalog number	Price
	in.	cm		
3	15½ x 26 x 14	39.4 x 66.0 x 35.6	GY-62510-00	
4	15½ x 37 x 14	39.4 x 94.0 x 35.6	GY-62510-10	
5	15½ x 48 x 14	39.4 x 122.0 x 35.6	GY-62510-20	

Jerricans with Handle, HDPE

Better stability with wider footprint

- Get accurate pouring with integral spout
- Don't lose the cap—attached via strap to keep in place
- Gallon and liter graduations



Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
6 (1½)	53B	GY-06069-10	1	
10 (2½)		GY-06069-20		
13 (3½)		GY-06069-25		
20 (5)		GY-06069-30		

[GY-63000-20](#) Replacement cap; PP, 53B. Pack of 12

UN-Approved Five-Gallon Container, HDPE

Easy dispensing where you need it

- UN- and FDA-compliant container
- Ship, store, or dispense samples easily

Liters (gallons)	Cap size	Catalog number	Qty/pk	Price/pk
20 (5)	70 mm	GY-06060-10	1	

[GY-06078-00](#) Plastic drum faucet, ¾" NPT(M), LDPE. Pack of 6



Container shown with plastic drum faucet

Carboys

Collapsible Containers, LDPE

Won't take up extra space—
easily expand when filled

- Stack neatly when empty, just unscrew cap and fill to expand
- Include polypropylene caps

Capacity L (gallon)	Cap size	Qty/ pk	Catalog number	Price/ pk
1 (1/4)	38 mm	12	GY-06100-10	
4 (1)			GY-06100-20	
10 (2 1/2)			GY-06100-30	
19 (5)		10	GY-06100-40	

**Dispensing Accessories for Carboys**

Get quicker sampling
with just a twist

Plastic Faucets, LDPE

Instant on/off control ensures fast, smooth flow. These faucets are manufactured of HDPE for greater strength and positive thread engagement; the stem is LDPE for 06078-00. Model 06079-00 has an EPDM gasket for a leakproof seal.



Fitting	Catalog number	Qty/pk	Price/pk
3/4" NPT(M)	GY-06078-00	5	
2" NPT(M)	GY-06079-00	1	

Spigots; PP, Tefzel® ETFE, and LDPE

Use these spigots with Nalgene carboys. Quick-action spigot 06070-22 includes an adapter for 1/4" to 3/8" ID tubing. A PTFE O-ring seated in the needle-type tank spigots ensures a leakproof seal.



Spigot OD	Material	Fitting	Catalog number	Qty/ pk	Price/ pk
Quick-action spigots; autoclavable					
5/8"	PP	1 1/8"-12 UNF	GY-06070-22	1	
	Tefzel ETFE	1 1/8"-12 UNF	GY-06070-23	1	
Needle-type tank spigots					
5/8"	LDPE	1/2" NPT(M)	GY-06070-20	1	
	LDPE	3/4" NPT(M)	GY-06070-30	1	

Closures for Bottles and Carboys
Fast and easy sampling

Flexible Top Works closures are available in solid (no ports), 2-port, and 3-port versions. Closures are polypropylene (PP). Platinum-cured silicone tubing is fused to a silicone stopper insert that forms a one-piece closure system. **Autoclavable** and leakproof.



Closure size	No. of insert ports	Port sizes	Recommended for	Catalog number	Price
38-430	Three	1/8", 1/8", 1/4"	Bottles/carboys	GY-13070-02	
53B	None (solid)	—	Bottles/jerricans	GY-13070-04	
53B	Two	1/4", 1/4"	Bottles/jerricans	GY-13070-06	
53B	Three	1/8", 1/4", 1/4"	Bottles/jerricans	GY-13070-08	
83B	None (solid)	—	Carboys	GY-13070-10	
83B	Two	1/4", 1/4"	Carboys	GY-13070-12	
83B	Three	1/4", 1/4", 3/8"	Carboys	GY-13070-14	

Wait!

**There's More at
ColeParmer.com**

Find more than
100,000 products
to choose from!