

Stainless Steel General-Purpose Carts

Move more items safely

- Optional guard rails keep items securely on the shelves
- All carts have swivel casters except 00946-12 which has two fixed and two swivel casters

Carts 00946-07 and -09 have 22-gauge steel shelves and 20-gauge steel legs; all other carts have 20-gauge steel shelves and 16-gauge steel legs.



00946-07

Shelf area (L x W)	No. of shelves	Capacity lb (kg)	Dimensions (L x W x H)	Catalog number	Price
General-purpose carts with guard rails					
24" x 15½" (61.0 x 39.4 cm)	3	200 (91)	27¾" x 16¾" x 33¼" (69.5 x 42.2 x 84.4 cm)	GY-00946-07	
27" x 18" (68.5 x 45.7 cm)	3	400 (181)	31¼" x 19¼" x 33¾" (79.3 x 48.9 x 85.7 cm)	GY-00946-08	
General-purpose carts without guard rails					
24" x 15½" (61.0 x 39.4 cm)	3	200 (91)	27¾" x 16¼" x 31⅞" (70.5 x 41.3 x 80.9 cm)	GY-00946-09	
27" x 18" (68.5 x 45.7 cm)	3	400 (181)	31¼" x 19" x 33¼" (79.3 x 48.3 x 84.4 cm)	GY-00946-10	
35" x 21" (88.9 x 53.3 cm)	3	400 (181)	39" x 22½" x 37½" (99 x 57.2 x 95.2 cm)	GY-00946-11	
49" x 21" (124.4 x 53.3 cm)	3	400 (181)	54" x 22½" x 37½" (137 x 57.2 x 95.2 cm)	GY-00946-12	

Stainless Steel Equipment Carts

No more carrying heavy equipment

- Optional drawer for storage of accessories
- Slide equipment on and off easily—lip on three sides to keep equipment in place

Models 00946-01 and -03 have a 22-gauge steel top and 20-gauge steel legs; models 00946-04 and -06 have 20-gauge steel tops and 18-gauge steel legs; models 00946-02 and 00946-05 have an 18-gauge steel top and 16-gauge steel legs.



00946-03

Shelf area (L x W)	No. of shelves	Capacity lb (kg)	Dimensions (L x W x H)	Catalog number	Price
Equipment carts					
18¾" x 16¾" (46.6 x 42.2 cm)	3	200 (91)	18¾" x 16¾" x 30" (47.6 x 42.5 x 76.2 cm)	GY-00946-01	
18¾" x 16¾" (46.6 x 42.2 cm)	3	400 (181)	18¾" x 16¾" x 30⅞" (47.6 x 42.5 x 78.4 cm)	GY-00946-02	
Equipment carts with drawer					
18¾" x 16¾" (46.6 x 42.2 cm)	3	200 (91)	18¾" x 16¾" x 30" (47.6 x 42.5 x 76.2 cm)	GY-00946-03	
18¾" x 16¾" (46.6 x 42.2 cm)	2	200 (91)	18¾" x 16¾" x 35" (47.6 x 42.5 x 88.9 cm)	GY-00946-04	
18¾" x 16¾" (46.6 x 42.2 cm)	3	400 (181)	18¾" x 16¾" x 30⅞" (47.6 x 42.5 x 78.4 cm)	GY-00946-05	
36" x 20" (91.4 x 50.8 cm)	3	400 (181)	36¾" x 20½" x 35" (92.4 x 51.1 x 88.9 cm)	GY-00946-06	

Heavy-Duty Stainless Steel Carts with Guard Rails

Don't worry about items slipping off shelves

- Keep loads on shelves with guard rails
- No damaged surroundings—bumpers on legs protect from damage

Shelves are made of 18-gauge steel; 1" (2.5 cm) square tube legs are made of ⅜" (3.2 mm)-thick steel.



00946-14

Shelf area (L x W)	No. of shelves	Capacity lb (kg)	Dimensions (L x W x H)	Catalog number	Price
27" x 18" (68.5 x 45.7 cm)	2	650 (295)	33¼" x 19½" x 34½" (84.4 x 49.5 x 87.6 cm)	GY-00946-13	
21" x 33" (53.3 x 83.8 cm)	2		38¾" x 22¾" x 37¼" (98.1 x 56.8 x 94.6 cm)	GY-00946-14	
	3			GY-00946-15	

Ergonomic Stainless Steel Carts

Reduces stress on arms and wrists

- Keep hands in comfortable position with vertical push handle
- Get a better grip with padded handles on 00946-18 and -19
- Minimize potential back strain with higher bottom shelf—20" from floor for less bending



00946-18

Shelf area (L x W)	No. of shelves	Capacity lb (kg)	Dimensions (L x W x H)	Catalog number	Price
24" x 15½" (61 x 39.4 cm)	2	300 (136)	28" x 16½" x 42⅞" (71.1 x 41.9 x 108.9 cm)	GY-00946-16	
27" x 18" (68.5 x 45.7 cm)	2		31⅞" x 19" x 42⅞" (79 x 48.3 x 108.9 cm)	GY-00946-17	
27" x 18" (68.5 x 45.7 cm)	2	500 (227)	31⅞" x 19" x 44½" (79 x 48.3 x 113 cm)	GY-00946-18	
35" x 21" (88.9 x 53.3 cm)	2		39⅞" x 22" x 44½" (99.4 x 55.9 x 113 cm)	GY-00946-19	

Carts

Specialty Laboratory Carts

No more wimpy carts—lifetime warranty guarantees long life

Watertight Pan Cart

- Contain spills with 18-gauge steel pans and 2" (5.1 cm) H sidewalls
- Prevent items from sliding with removable rubber mats



Capacity lb (kg)	Dimensions (L x W x H)		Catalog number	Price
	in.	cm		
400 (181)	34 $\frac{7}{8}$ x 19 x 36 $\frac{3}{4}$	89 x 48 x 92	GY-09389-00	

Flexi-Bin Carts

- Secure and protect fragile items for transport in removable, corrosive-resistant polyethylene bin
- Durable epoxy-coated steel frame, molded plastic handle grips, and two locking casters



Capacity lb (kg)	Dimensions (L x W x H)		Catalog number	Price
	in.	cm		
400 (181)	34 $\frac{7}{8}$ x 19 x 36 $\frac{3}{4}$	89 x 48 x 92	GY-09380-02	

[GY-09380-04](#) Replacement bin,
22 $\frac{1}{4}$ "W x 17"H x 10"D (57 x 43 x 25 cm)

Chemical Cart

- Withstand occasional spills of acids and caustics
- Laminate shelves resist chemicals, heat, and stains



Capacity lb (kg)	Dimensions (L x W x H)		Catalog number	Price
	in.	cm		
400 (181)	34 $\frac{7}{8}$ x 19 x 36 $\frac{3}{4}$	89 x 48 x 92	GY-09374-00	

Glassware Cart

- Ideal for washing and drying labware—drip pan catches liquid and has plastic-coated shelves
- Customize with large or small replacement baskets



No. of baskets held	Dimensions (L x W x H)		Catalog number	Price
	in.	cm		
2 large	34 $\frac{7}{8}$ x 19 x 36 $\frac{3}{4}$	89 x 48 x 92	GY-09378-00	
4 small			GY-09382-00	

[GY-09380-00](#) Large replacement basket,
32"L x 18"W x 4"H (81 x 46 x 10 cm)

[GY-09384-00](#) Small replacement basket,
16"L x 18"W x 4"H (41 x 46 x 10 cm)

Carboy Caddy

No more stressful moving of heavy containers

- Tip cart, slide carboy on cart, and push—it's stress-free moving
- Get easy navigation with 360° swivel casters
- Constructed of phenolic base with melamine laminate surface



Capacity lb (kg)	Platform dimensions		Handle height	Catalog number	Price
	in.	cm			
360 (163)	20 x 23 $\frac{5}{8}$	51 x 60	38 $\frac{1}{2}$ " (98 cm)	GY-09378-05	

Utility Carts

No more bacterial, mold, mildew, or fungal growth

- Antimicrobial treatment protects shelves
- No corrosion or rusting—made of strong polymer
- Unique open-grid shelf design
- Easily removable shelves allow for cleaning



Capacity lb (kg)	No. of shelves	Shelf area (L x W)	Overall height		Catalog number	Price
			in.	cm		
600 (272)	2	36" x 24" (91 x 61 cm)	40	102	GY-47300-22	
600 (272)	3	30" x 18" (76 x 46 cm)	40	102	GY-47300-24	
		36" x 24" (91 x 61 cm)			GY-47300-26	

Space-Saving Folding Carts

Won't take up extra space—
easily folds up

- Fold and unfold in one easy motion!
- Takes up minimal lab space
- Durable aluminum and polyethylene construction
- Swivel casters for easy movement, two lockable casters for securing

Folds
without
tools!



Capacity lb (kg)	No. of shelves	Dimensions (L x W x H)		Catalog number	Price
		in.	cm		
350 (159)	3	30 x 15 x 35½, folds to 9"W	76 x 38 x 90, folds to 23 cm W	GY-47600-00	
		40 x 20½ x 35½, folds to 10"W	102 x 52 x 90, folds to 25 cm W	GY-47600-10	

Low-Cost Service Carts

Don't worry about cleanliness—
steam-cleanable design

- Easily contain spills with seamless, 2¾" (7 cm) deep shelves
- Molded, higher handle for less strain on arms and back
- Large rubber casters provide easy mobility—two casters swivel



Capacity lb (kg)	No. of shelves	Shelf area (L x W)	Dimensions (L x W x H)		Catalog number	Price
			in.	cm		
400 (181)	2	30½" x 16½" (77 x 41 cm)	34½" x 17¾" x 32¾" (87 x 44 x 83 cm)		GY-09348-00	
		36¼" x 24¼" (92 x 62 cm)	40½" x 25¾" x 32¾" (103 x 65 x 83 cm)		GY-09348-10	

Lockable Specialty Carts

No more stolen supplies—
lockable cabinets for safe storage

- Get two uses in one—convenient height for use as mobile workbench

X-tra™ Instrument Cart

- Keep supplies at hand and minimize clutter with pull-out drawer
- Easy-to-clean shelves and cabinet
- Swivel casters for smooth transport



Capacity lb (kg)	Dimensions (L x W x H)		Catalog number	Price
	in.	cm		
300 (136)	40½ x 20 x 37½	103 x 51 x 96	GY-47205-30	

X-tra Equipment Cart

- Keep items concealed with three enclosed sides and one lockable cabinet
- Handles up to 100 lb (45.4 kg) per shelf
- Swivel casters for easy movement



Capacity lb (kg)	Dimensions (L x W x H)		Catalog number	Price
	in.	cm		
300 (136)	40½ x 20¾ x 37½	103 x 53 x 96	GY-47205-40	

Trades Cart

- Convenient storage areas to hold a variety of items
- Mounting bracket for electrical cord or power strip, storage bin for small parts, and four compartmentalized boxes
- Two rigid casters and two swivel casters with locks



Capacity lb (kg)	Dimensions (L x W x H)		Catalog number	Price
	in.	cm		
350 (160)	38 x 19¼ x 33½	97 x 49 x 85	GY-47205-50	

Carts

Flat Shelf Carts

No more rusting, denting, or chipping

- Easily slide heavy items on and off flat top shelf
- Large, smooth-gliding casters won't mar hard surfaces
- Foam plastic construction for strength, quiet operation, and easy cleaning



Capacity lb (kg)	No. of shelves	Shelf area (L x W)	Dimensions (L x W x H)	Catalog number	Price
400 (182)	2	28" x 15¾" (71 x 40 cm)	38" x 19" x 33½" (97 x 48 x 85 cm)	GY-09348-30	
		34" x 22¼" (86 x 57 cm)	44" x 25½" x 33½" (112 x 65 x 85 cm)	GY-09348-40	

Utility Carts

No more muscle strain

- Ergonomic handles and easy-gliding casters for smooth handling
- Constructed of rugged, webbed-polypropylene shelves and aluminum supports
- Model 47205-20 has side panels on three sides to conceal contents
- 300 lb (135 kg) capacity



Description	No. of shelves	Dimensions (L x W x H)		Catalog number	Price
		in.	cm		
Cart, open sides	3	40½ x 20 x 37½	103 x 51 x 96	GY-47205-00	
Cart, enclosed on three sides				GY-47205-20	

General-Purpose Reversible Tray Carts

Versatile—like two carts in one

- Top tray inverts for either flat-top use or shelf with sides
- Keep items from sliding off shelves with deep 4" (10.2 cm) trays
- No snagged clothing—carts have no exposed nuts, bolts, or rivets
- Enamel-coated steel with four casters, two swivel



Capacity lb (kg)	No. of shelves	Shelf area (L x W)	Overall height		Catalog number	Price
			in.	cm		
500 (227)	2	30" x 16" (76 x 41 cm)	36	91	GY-09367-20	
		36" x 24" (91 x 61 cm)			GY-09365-20	

Welded Service Carts

Holds just about anything—1000-lb (424-kg) capacity

- Choose from flat nonlipped or lipped shelving option—lipped shelves keep contents from rolling off
- Constructed of 12-gauge heavy-duty steel
- Two swivel casters provide easy maneuvering
- 1000 lb (424 kg) capacity



No. of shelves	Shelf area (L x W)	Overall height		Flat shelves	Lipped shelves	Price
		in.	cm	Cat. no.	Cat. no.	
2	18" x 32" (46 x 81 cm)	35	89	GY-47355-74	GY-47355-82	
	24" x 36" (61 x 91 cm)			GY-47355-76	GY-47355-84	
3	18" x 32" (46 x 81 cm)	35	89	GY-47355-78	GY-47355-86	
	24" x 36" (61 x 91 cm)			GY-47355-80	GY-47355-88	

Wait!

**There's More at
ColeParmer.com**

Browse **100,000+** products plus find
Technical Resources to help you choose

Lab Dolly

Protect against bacteria, mold, and mildew

- Easily transport heavy loads throughout your lab or facility
- Features a polymer platform with stainless steel handle and frame for corrosion resistance and long life
- Recessed-tray design contains spills and wipes clean easily



Capacity lb (kg)	Platform dimensions (L x W)	Dimensions (L x W x H)	Catalog number	Price
150 (68)	22" x 33" (56 x 84 cm)	22" x 33" x 41" (56 x 84 x 104 cm)	GY-09378-01	

Deep-Ledge Utility Carts

No more spills

- Tough, high-density polyethylene provides excellent chemical resistance
- Two- and three-shelf carts manufactured with Microban®, an antimicrobial product that protects from bacteria, mold, and mildew
- 400-lb (181 kg) capacity



47300-18

Shelves are smooth and easily wipe clean. The middle shelf of the three-shelf model adjusts at 1" (2.5 cm) increments. All models come with 5"-dia (12.7 cm) swivel casters.

No. of shelves	Shelf area (L x W)	Shelf clearance	Height in.	Height cm	Color	Catalog number	Price
2	32¾" x 21½" (83 x 55 cm)	38" (97 cm)	41	104	Gray	GY-47300-16	
3		16" (41 cm)			Blue	GY-47300-18	

Heavy-Duty Service Carts with Optional Bins

Shelves adjust for different load sizes

- Impact-resistant polyolefin shelves supported by 4"-dia (10.2 cm) rubber casters
- Three-shelf models include adjustable middle shelf
- Shipped unassembled



Cart 47300-10 shown with optional wastebasket 47300-50



Capacity lb (kg)	No. of shelves	Shelf area (L x W)	Dimensions (L x W x H)	Catalog number	Price
300 (136)	2	28½" x 18½" (73 x 47 cm)	33¾" x 21½" x 33¼" (86 x 55 x 84 cm)	GY-47300-00	
400 (181)	3			GY-47300-10	
300 (136)	2	34½" x 24½" (88 x 62 cm)	39½" x 27" x 33¼" (100 x 69 x 84 cm)	GY-47300-02	
400 (181)	3			GY-47300-12	

GY-47300-50 Wastebasket (includes holder) for carts 47300-00 and -10 only. Measures 18½"W x 23¼"H x 12"D (47 x 59 x 31 cm)

GY-47300-60 Utility bin (includes holder) for carts 47300-00 and -10 only. Measures 18½"W x 10½"H x 12¾"D (47 x 27 x 32 cm)

Nonconductive Instrument Carts

No more static discharge

- Shelves are nonconductive and made of MDF with a nonslip surface

Carts 47352-00 and -20 are cushion-loaded and ideally suited for transporting delicate equipment or sensitive devices.

Carts 47355-24 and -26 feature a UL-listed, 15-amp power strip with six grounded outlets, a 15-ft (4.6-m) power cord, circuit breaker, and switch. Generous storage drawer with cylinder lock provides added security.



Cart 47355-24 with outlet strip and drawer

Specifications

Capacity: 600 lb (272 kg) **Overall height:** 34" (86 cm)



47355-24, -26 only

Description	No. of shelves	Shelf area (L x W)	Shelf clearance	Catalog number	Price
Basic cart	2	Top: 36" x 22" (91 x 56 cm) Bottom: 30" x 20" (76 x 51 cm)	12" (31 cm)	GY-47352-00	
			16" (41 cm)	GY-47352-20	
Cart with outlet strip and drawer			12" (31 cm)	GY-47355-24	
			16" (41 cm)	GY-47355-26	

Heavy-Duty Two-Shelf Wire Utility Carts

Transport up to 800 lb (363 kg)!

Carts have four 5"-dia (12.7 cm) rubber swivel casters.

Capacity lb (kg)	No. of shelves	Shelf size (L x W)	Overall height	Catalog number	Price
800 (363)	2	48" x 18" (122 x 46 cm)	39" (99 cm)	GY-47521-10	
		36" x 24" (91 x 61 cm)		GY-47521-30	
		48" x 24" (122 x 61 cm)		GY-47521-40	
		60" x 24" (152 x 61 cm)		GY-47521-50	

