

Containers

Cole-Parmer® Disposable Lab Containers, PPCO

No more wasted space—containers stack easily

- Durable for longer use—stand up to boiling, freezing, and most chemicals
- Know what's inside—translucent containers come with tight-fitting lids



Cole-Parmer®

Capacity mL (oz)	Catalog number	Qty/pk	Price/pk
120 (4)	GY-06727-00	250	
240 (8)	GY-06727-10	250	
480 (16)	GY-06727-20	100	
960 (32)	GY-06727-30	100	
2000 (64)	GY-06727-35	50	
2500 (83)	GY-06727-42	25	
5000 (170)	GY-06727-50	25	

Utility Boxes, Styrene-Acrylonitrile

Stack to save space

- Store dry materials easily while knowing what's inside
- Freezable
- Easy to clean
- Include snap-on lids for tight storage



Capacity mL (oz)	Dimensions (L x W x H)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
325 (10)	5 1/16 x 2 15/16 x 2 5/16	12.9 x 7.5 x 5.9	GY-63524-20	4	
500 (16)	7 1/2 x 6 1/8 x 1 7/16	19.1 x 15.6 x 3.7	GY-63524-30	2	
1000 (32)	5 x 2 7/8 x 6 7/8	12.7 x 7.3 x 17.5	GY-63524-50		
2000 (64)	7 1/2 x 6 1/8 x 3 11/16	19.1 x 15.4 x 9.4	GY-63524-60	1	

Containers, HDPE

Be confident lid is tight

- Optional snap-tight lids feature positive lock for tight seal
- Can be frozen or chemically sterilized
- Easy carrying grip handles
- NSF certified; FDA-compliant



Capacity L (qt)	Dimensions (H x OD)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
2 (2)	4½ x 7	11.4 x 17.8	GY-63530-00	12	
4 (4)	7½ x 7	19.1 x 17.8	GY-63530-02		
6 (6)	8 x 8¾	20.3 x 22.2	GY-63530-04		
8 (8)	10 x 8¾	25.4 x 22.2	GY-63530-06	12	
12 (12)	12 x 10	30.5 x 25.4	GY-63530-08		
16 (16)	9½ x 12	24.1 x 30.5	GY-63530-10		
20 (20)	13½ x 12	34.3 x 30.5	GY-63530-12	6	
24 (24)	15 x 12	38.1 x 30.5	GY-63530-14		
32 (32)	12 x 15½	30.5 x 39.4	GY-63530-16		
40 (40)	15½ x 15½	39.4 x 39.4	GY-63530-18		
Snap-tight, positive-lock lids					
For 2- and 4-L containers			GY-63530-50	12	
For 6- and 8-L containers			GY-63530-51		
For 12-L container			GY-63530-52		
For 16-, 20-, 24-L containers			GY-63530-53	6	
For 32- and 40-L containers			GY-63530-54		

Round Containers with Handles, PC

Leakproof design for easy transport

- Easily identify contents with translucent design
- Know volume with quart and liter graduations
- Withstand temperatures from -40 to 212°F (-40 to 100°C)
- NSF certified



Capacity L (qt)	Dimensions (H x overall dia)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
1 (1)	5 x 6 3/16	12.7 x 15.7	GY-63571-00	12	
2 (2)	4 3/4 x 8 3/16	10.8 x 20.8	GY-63571-02	12	
4 (4)	8 7/16 x 8 3/16	21.4 x 20.8	GY-63571-04	12	
6 (6)	7 3/4 x 10	19.7 x 25.4	GY-63571-06	12	
8 (8)	10 3/4 x 10	27.3 x 25.4	GY-63571-08	12	
12 (12)	8 1/4 x 13 3/8	21.0 x 35.2	GY-63571-10	6	
18 (18)	11 7/8 x 13 3/8	30.2 x 35.2	GY-63571-12	6	
22 (22)	15 x 13 3/4	38.1 x 34.9	GY-63571-14	6	
Snap-on lids					
For 1-L containers, translucent			GY-63571-20	12	
For 2- and 4-L containers, green			GY-63571-22	12	
For 6- and 8-L containers, red			GY-63571-24	12	
For 12-, 18-, and 22-L containers, red			GY-63571-26	6	

Storage Containers with Handles, HDPE and PP

No more difficult lifting

- Easy-lifting grip handles, 22-qt containers have removable handles
- Keep track of contents with “date-content” label area
- Feature graduations in quarts and liters and optional tight fitting lids
- Resist temperatures from –40 to 150°F (–40 to 66°C)
- Dishwasher safe, NSF certified



HDPE
63511-01



PP
63512-13

Lids sold separately



Capacity L (qt)	Dimensions (H x OD)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
White high-density polyethylene (HDPE) containers					
2 (2)	4 x 8½	10.2 x 21.6	GY-63511-01	12	
4 (4)	7¾ x 8½	19.7 x 21.6	GY-63511-03		
6 (6)	7⅝ x 10	19.4 x 25.4	GY-63511-05	12	
8 (8)	10⅝ x 10	27.0 x 25.4	GY-63511-07		
12 (12)	8⅝ x 13⅝	20.6 x 33.3	GY-63511-09	6	
18 (18)	11⅞ x 13⅝	30.2 x 33.3	GY-63511-11		
22 (22)	14 x 13⅝	35.6 x 33.3	GY-63511-13		
Clear polypropylene (PP) containers					
2 (2)	4 x 8½	10.2 x 21.6	GY-63512-01	12	
4 (4)	7¾ x 8½	19.7 x 21.6	GY-63512-03		
6 (6)	7⅝ x 10	19.4 x 25.4	GY-63512-05	12	
8 (8)	10⅝ x 10	27.0 x 25.4	GY-63512-07		
12 (12)	8⅝ x 13⅝	20.6 x 33.3	GY-63512-09	6	
18 (18)	11⅞ x 13⅝	30.2 x 33.3	GY-63512-11		
22 (22)	14 x 13⅝	35.6 x 33.3	GY-63512-13		
Yellow low-density polyethylene (LDPE) lids					
For 2- and 4-L containers			GY-63511-53	12	
For 6- and 8-L containers			GY-63511-54	12	
For 12- 18- 22-L containers			GY-63511-55	6	

DON'T RISK IT!

Our chemical compatibility database provides important information on chemical resistance—check it out at ColeParmer.com/ChemChart



Keep track of contents by recording date and material on lid.

Round Dipping Baskets, PP

Easily remove contents from solution

- ⅛" perforations for liquid release
- Easy-grip handle for ease of handling
- Use with crocks sold at right
- NSF certified



Round dipping basket
06717-50
(shown with heavy crock)

Heavy Crocks, HDPE

Replace heavier, older style crocks

- Sturdy design—one piece molded HDPE
- Use with dipping baskets sold at left for easy item retrieval
- Include cover for easy storage
- NSF certified; USDA approved



Heavy crock
06724-60
(cover included)
shown with
round dipping
basket

Fits crock model no.	Dimensions (H x ID)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
06724-20	8 x 8	20.3 x 20.3	GY-06717-20	1	
06724-40	10 x 9	25.4 x 22.9	GY-06717-30		
06724-60	13 x 9	33.0 x 22.9	GY-06717-50		
06724-80	17 x 11½	43.2 x 29.2	GY-06717-70		

Capacity L (gal.)	Dimensions (H x OD)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
12 (3.25)	11¼ x 10½	29.8 x 25.7	GY-06724-20	1	
15 (4)	13¼ x 12	33.7 x 30.5	GY-06724-40		
22 (5.75)	16½ x 11	41.9 x 27.9	GY-06724-60		
40 (10.25)	21 x 13	53.3 x 33.0	GY-06724-80		

Square Utility Pan, 304 SS (medium gauge)

Holds up to tough demands

- Sturdy, stainless steel construction for longer life
- NSF certified

Capacity L (qt)	Dimensions (L x W x H)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
1.25 (1¼)	6½ x 6½ x 2½	15.6 x 15.6 x 7.3	GY-07286-01	1	



Containers

Utility Basins, HDPE

Safely transport what you need

- Molded-in handles for easy carrying
- Won't take up extra space—tapered design for nesting when not in use



Capacity L (qt)	Outside dimensions (L x W x H)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
10 (10)	14½ x 12½ x 5½	35.9 x 30.8 x 14.3	GY-06011-15	1	
15 (15)	21½ x 17½ x 5½	54.3 x 43.5 x 13.0	GY-06011-30	1	

Boxes and Lids, HDPE

Not your typical box

- Multipurpose box and system configure to your needs
- Add colander for easy straining of contents
- Optional lid doubles as tray

Box 06017-16 shown with optional colander 06017-50 and optional lid 06017-22



Capacity L (gal.)	Outside dimensions (L x W x H)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
Boxes					
8 (2)	18 x 12 x 3½	45.7 x 30.5 x 8.9	GY-06017-10	1	
13 (3½)	18 x 12 x 6	45.7 x 30.5 x 15.2	GY-06017-12		
19 (5)	18 x 12 x 9	45.7 x 30.5 x 22.9	GY-06017-14		
26.5 (7)	18 x 26 x 3½	45.7 x 66 x 8.9	GY-06017-15		
32 (8½)	18 x 26 x 6	45.7 x 66 x 15.2	GY-06017-16		
47 (12½)	18 x 26 x 9	45.7 x 66 x 22.9	GY-06017-18		
92 (24½)	18 x 26 x 12	45.7 x 66 x 30.5	GY-06017-19		
81 (21½)	18 x 26 x 15	45.7 x 66 x 38.1	GY-06017-20		
Lids					
For 06017-10, -12, -14			GY-06017-21	1	
For 06017-15 through -20			GY-06017-22		
Colander					
For 06017-16 -18 -19 -20			GY-06017-50	1	

Basins and Cover, HDPE

Get organized—stackable design

- Keep contents covered with optional lid
- Nest for easy storage when not in use
- Stand up to temperatures from –40 to 180°F (–40 to 83°C)
- Contain no BPA's; NSF certified

06019-10 shown with cover 06019-01



Outside dimensions (L x W x H)		Catalog number	Qty/ pk	Price/pk
in.	cm			
Basins				
21 x 15½ x 5	53.3 x 39.4 x 12.7	GY-06019-10	6	
21 x 15½ x 7	53.3 x 39.4 x 17.8	GY-06019-20	3	
Cover				
For basins 06019-10, -20		GY-06019-01	6	

Collapsible Basket, PP

Won't take up extra space

- Fold down when not in use to just 2¾" high
- Easy carrying with molded in handles
- Stainless steel pin hinges stand up to your demands



06818-70 with optional casters

Collapsed position

Basket color	Capacity		Dimensions (L x W x H)		Catalog number	Qty/ pk	Price/pk
	cu ft	lb	in.	cm			
Blue	1.5	70	20¾ x 13¾ x 10¾	52.7 x 34.9 x 27.3	GY-06818-70	1	
Red					GY-06818-72		

Graduated Sample Containers, PP

Leakproof for better sampling

- Easier sampling—graduated for measuring directly
- Durable design—translucent, unbreakable, and chemical resistant
- Sterile containers are individually wrapped



Capacity mL (oz)	Description	Catalog number	Qty/ pk	Price/pk
118 (4)	Sterile, no tab	GY-06049-40	100	
118 (4)	Nonsterile with white tab	GY-06049-50	500	

Rectangular Mobile Tubs, LDPE

No more strained backs

- Roll even heavy loads around with ease on rugged caster frame and wheels
- Easy clean-up—remove tub and wheels for washing
- Include polyethylene cover
- NSF certified and USDA approved



Capacity L (gal.)	Dimensions (L x W x H*)		Catalog number	Price/ ea
	in.	cm		
60 (15)	32 $\frac{3}{8}$ x 15 $\frac{3}{8}$ x 12 $\frac{1}{4}$	82.2 x 39.1 x 31.1	GY-06745-20	
106 (28)	26 $\frac{3}{4}$ x 20 $\frac{1}{2}$ x 16 $\frac{3}{4}$	67.9 x 52.1 x 42.5	GY-06745-40	
204 (54)	35 $\frac{1}{2}$ x 26 x 20 $\frac{1}{4}$	90.2 x 66.0 x 51.4	GY-06745-60	

Accessories

[GY-60104-02](#) HDPE bag, 24" x 33" (61 x 84 cm)

[GY-60104-04](#) HDPE bag, 30" x 37" (76 x 94 cm)

[GY-60104-06](#) HDPE bag, 33" x 40" (84 x 102 cm)

Cole-Parmer® Sample Containers, PP

Resist most chemicals

- Chemical resistant to acids, bases, and most organic solvents
- Easy read graduations in mL and oz
- Lids fit tightly, sold separately



Capacity mL (oz)	Dimensions (H x OD)		Catalog number	Qty/ pk	Price/ pk
	in.	cm			
Sample containers					
120 (4)	2 ⁹ / ₁₆ x 2 ⁵ / ₈	6.5 x 6.7	GY-06100-12	500	
240 (8)	4 ¹ / ₁₆ x 2 ⁵ / ₈	10.3 x 6.7	GY-06100-14	500	
Lids					
—	⁵ / ₁₆ x 2 ³ / ₄	0.8 x 7.0	GY-06100-16	500	

Heavy-Duty, Round Mobile Containers, HDPE

Don't buy a separate cart

- Handle liquids and solids with ease
- HDPE cover keeps contents intact during transport
- Container and casters are removable for easier cleaning
- NSF certified and USDA approved



Capacity L (gal.)	Dimensions (H x OD)		Catalog number	Qty	Price/ ea
	in.	cm			
100 (26)	29 $\frac{3}{4}$ x 20 $\frac{3}{4}$	75.6 x 52.7	GY-06740-20	1	
197 (52)	39 x 22	99.1 x 55.9	GY-06740-30	1	

Containers

Hinged-Lid Sample Containers, PP

Stop losing lids

- Attached lid is always close at hand
- Won't leak or spill when closed
- FDA-compliant
- Autoclavable



Capacity mL (oz)	Dimensions (H x OD)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
7.5 (¼)	5/8 x 1½	1.6 x 3.8	GY-66098-00	2500	
15 (½)	¾ x 1½	1.9 x 3.8	GY-66098-10	2000	
22.5 (¾)	1 x 1½	2.5 x 3.8	GY-66098-20	1500	
30 (1)	1½ x 1½	3.8 x 3.8	GY-66098-30	1300	
60 (2)	1½ x 2	4.1 x 5.1	GY-66098-40	750	
120 (4)	3¼ x 2	8.1 x 5.1	GY-66098-50	300	

Accessories

[GY-06098-01](#) Labels for containers 66098-00 through -30.
Pack of 1000

[GY-06098-02](#) Labels for containers 66098-40 and -50. Pack of 1000

Wide-Mouth Sample Containers, PP

Mix and store in one vessel

- Easily attach to blenders for in-container mixing
- Models 06101-31 and -40 fit blenders with Mason jar threads
- Durable for general storage—chemical-resistant polypropylene
- Autoclavable



Capacity mL (oz)	Cap size	Catalog number	Qty/ pk	Price/ pk
30 (1)	43 mm	GY-06101-11	12	
60 (2)	53 mm	GY-06101-21		
120 (4)	70 mm	GY-06101-31		
240 (8)	89 mm	GY-06101-40	6	
480 (16)	89 mm	GY-06101-50		
960 (32)	120 mm	GY-06101-60		
1190 (40)	120 mm	GY-06101-70		

Snap-Seal™ Sample Containers, PP

Won't spill no matter the orientation

- Lids snap shut for leakproof fit
- Locking lids eliminate inadvertent opening
- Autoclavable



Unique seal won't spill or leak.



62261-30

62261-00

Capacity mL (oz)	Dimensions (H x OD)		Catalog number	Qty/ pk	Price/pk
	in.	cm			
13 (⅞)	3¼ x 5/8	9.5 x 1.6	GY-62261-00	500	
45 (1½)	3¾ x 1⅝	8.6 x 2.9	GY-62261-10	400	
120 (4)	3½ x 1¾	8.9 x 4.4	GY-62261-20	200	
240 (8)	3 x 3⅝	7.6 x 7.9	GY-62261-30	100	
300 (10)	4¾ x 2½	11.1 x 6.4	GY-62261-40	100	

Straight-Sided Sample Jars, Glass

Don't strain to get samples

- Get easy access to contents with wide-mouth design
- Store light-sensitive items safely in amber jars
- Lids are included
- Made of Type III soda lime glass



Capacity mL (oz)	Cap size	Clear jars			Amber jars		
		Cat. no.	Qty/cs	Price/cs	Cat. no.	Qty/cs	Price/cs
60 (2)	53-400	GY-34700-01	24		GY-34700-51	24	
120 (4)	53-400	GY-34700-03			GY-34700-53		
180 (6)	63-400	GY-34700-05			—		
240 (8)	70-400	GY-34700-07	12		GY-34700-57	12	
480 (16)	89-400	GY-34700-09			GY-34700-59		
960 (32)	89-400	GY-34700-11			GY-34700-61		

Jars, PFA

Samples won't stick to surface

- Nonstick surface is easy to clean
- Resist corrosive liquids, solvents, acids, sensitive reagents, and staining solutions
- Use 2200-mL size for EPA TCLP soil sampling (40 CFR, Part 261 et. atl.)
- Autoclavable



Capacity mL (oz)	Dimensions (H x OD)		Neck size	Catalog number	Price/ea
	in.	cm			
60 (2)	1 1/8 x 2	4.7 x 5.1	51 mm	GY-34607-58	
120 (4)	2 x 2 1/2	5.1 x 6.4	66 mm	GY-34607-60	
240 (8)	4 1/4 x 3 1/4	10.8 x 8.3	84 mm	GY-34607-62	
360 (12)	3 1/2 x 3 1/4	8.8 x 8.3	84 mm	GY-34607-64	
480 (16)	3 3/4 x 2 1/2	9.4 x 6.4	65 mm	GY-34607-66	
1000 (32)	6 x 4 1/4	15.2 x 10.7	107 mm	GY-34607-68	
2000 (64)	10 1/2 x 4	26.9 x 10.4	105 mm	GY-34607-70	
2200 (70)	11 1/2 x 4 1/4	29.3 x 10.7	107 mm	GY-34607-72	

Wide-Mouth Sample Containers, PS

Won't break like glass



- Get glass-like clarity with more rigidity
- Included screw caps keep samples secure
- Containers are polystyrene (PS); lids are polypropylene (PP) with polyethylene foam lines

Capacity mL (oz)	Dimensions (H x OD)		Cap size	Catalog number	Qty/ pk	Price/pk
	in.	mm				
60 (2)	2 x 1 3/4	53 x 46	53 mm	GY-62263-20	48	
120 (4)	2 3/4 x 2	68 x 53	58 mm	GY-62263-30	36	
240 (8)	3 1/2 x 2 1/2	89 x 65	89 mm	GY-62263-40	36	
480 (16)	3 3/4 x 3 1/2	96 x 89	89 mm	GY-62263-50	24	
960 (32)	4 3/4 x 3 3/8	120 x 98	120 mm	GY-62263-60	24	

Cole-Parmer® Cylindrical Jars, HDPE

Get easy sample access

- Wide neck for easy powder sampling
- Include screw caps to keep samples intact
- Won't break—made of high-density polyethylene for sturdy storage



Capacity mL (oz)	Dimensions (H x OD)		Catalog number	Qty/ pk	Price/pk
	in.	mm			
70 (2.4)	2 5/16 x 2	60 x 50	GY-63100-10	10	
120 (4)	2 3/4 x 2 1/4	71 x 56	GY-63100-12		
250 (8)	3 3/4 x 2 3/4	94 x 69	GY-63100-14		
500 (16)	4 1/4 x 3 3/8	107 x 86	GY-63100-16		
1000 (32)	5 x 4 3/8	128 x 111	GY-63100-18	10	
1500 (51)	7 7/16 x 4 3/8	182 x 111	GY-63100-20		
2000 (64)	9 1/4 x 4 3/8	235 x 111	GY-63100-22		

Collapsible Containers, LDPE

Easily expand when filled to save space

- Stack neatly when empty, just unscrew cap and fill to expand
- Include polypropylene caps



Capacity L (gal.)	Cap size	Catalog number	Qty/ pk	Price/ pk
1 (1/4)	38 mm	GY-06100-10	12	
4 (1)		GY-06100-20		
10 (2 1/2)		GY-06100-30	10	
19 (5)		GY-06100-40		