

High-Form Crucibles and Covers, Porcelain

Effortless cleaning with glazed surfaces

High-Form Crucibles are glazed inside and out except for outside bottom and rim.

Covers are glazed inside and out except for rim. **Autoclavable.**

Capacity (mL)	Top outside dia (mm)	Height (mm)	Crucibles					
			Cat. no.	Qty/pk	Price/pk	Cat. no.	Qty/cs	Price/cs
1.3	15	15	GY-17803-00	6		GY-17803-20	24	
5	24	20	GY-17803-02	12		GY-17803-22	72	
10	31	26	GY-17803-04	12		GY-17803-24	72	
15	35	29	GY-17803-06	12		GY-17803-26	72	
30	43	37	GY-17803-08	12		GY-17803-28	36	
40	50	40	GY-17803-10	12		GY-17803-30	36	
50	53	43	GY-17803-12	12		GY-17803-32	24	
100	65	54	GY-17803-14	6		GY-17803-34	24	
250	88	72	GY-17803-16	6		GY-17803-36	36	

Fits crucible catalog number	Inside rim dia (mm)	Covers					
		Cat. no.	Qty/pk	Price/pk	Cat. no.	Qty/cs	Price/cs
17803-02 and -22	26	GY-17811-00	12		GY-17811-20	72	
17803-04 and -24	34	GY-17811-02			GY-17811-22		
17803-06 and -24	38	GY-17811-04			GY-17811-24		
17803-08 and -24	45	GY-17811-06			GY-17811-26		
17803-10 and -24	54	GY-17811-08	6		GY-17811-28	72	
17803-12 and -24	66	GY-17811-10			GY-17811-30	60	
17803-14 and -24	90	GY-17811-14			GY-17811-34	36	



Wide-Form Crucibles and Covers, Porcelain

Easier sample placement with wider mouth

Wide-Form Crucibles are glazed inside and out except for outside bottom and rim. **Covers** are glazed inside and out except for rim. **Autoclavable.** Case quantities available at ColeParmer.com

Capacity (mL)	Top outside dia (mm)	Height (mm)	Crucibles			Covers		
			Cat. no.	Qty/pk	Price/pk	Cat. no.	Qty/cs	Price/cs
8	32	20	GY-17804-02	12		GY-17811-02	12	
17	41	26	GY-17804-04	12		GY-17811-06	12	
30	50	31	GY-17804-06	12		GY-17811-08	6	
50	61	37	GY-17804-08	12		GY-17811-10	6	
100	76	46	GY-17804-10	12		GY-17811-12	6	
150	84	52	GY-17804-12	1		GY-17811-14	6	
250	102	60	GY-17804-14	1		GY-17811-16	6	



Cole-Parmer® Crucibles, Quartz

Won't crack under extreme conditions

- Get low expansion, excellent dimensional stability, and thermal shock resistance
- Withstand temperatures up to 1922°F (1050°C)
- Crucible and lid sold separately



Capacity (mL)	Top outside dia (mm)	Height (mm)	Crucibles		Covers	
			Catalog number	Price/ea	Catalog number	Price/ea
15	41	25	GY-34526-30		GY-34528-30	
25	47	28	GY-34526-32		GY-34528-32	
50	57	37	GY-34526-34		GY-34528-34	
80	67	45	GY-34526-36		GY-34528-36	
100	70	48	GY-34526-38		GY-34528-38	
150	82	51	GY-34526-40		GY-34528-40	

NEW



Crucibles

Crucible and Cover, Quartz

Stand up to the heat

- Withstands temperatures to 1922°F (1050°C)
- Ensure chemical purity for subsequent analysis



Description	Capacity (mL)	OD x H (mm)	Catalog number	Price/ea
Crucible and cover	50	49 x 50	GY-34526-50	
Crucible only			GY-34526-52	

[GY-34526-54](#) Replacement quartz crucible cover

High-Form Crucibles and Covers, High Alumina

Won't give into chemicals

- Highly resistant to chemicals, alkalines, and fluxes
- Stands up to the heat—up to 3182°F (1750°C)
- Inactive in hydrogen and carbonaceous atmosphere



Covers sold separately

Capacity (mL)	Top outside dia (mm)	Height (mm)	Crucibles		Covers	
			Catalog number	Price/ea	Catalog number	Price/ea
10	29	31	GY-17804-01		GY-17804-11	
20	38	35	GY-17804-03		GY-17804-13	
50	48	47	GY-17804-05		GY-17804-15	
100	58	61	GY-17804-07		GY-17804-17	
250	76	93	GY-17804-09		GY-17804-19	

Cole-Parmer® High-Form Crucibles and Lids, Nickel

Get more strength—seamless

- Constructed from one piece of 24 gauge (0.025") nickel
- Polished inside and out
- Withstand temperatures to 1472°F (800°C)



Covers sold separately



Capacity (mL)	Top outside dia (mm)	Bottom outside dia (mm)	Height (mm)	Crucibles		Covers	
				Catalog number	Price/ea	Catalog number	Price/ea
15	31	26	31	GY-07281-50		GY-07281-60	
25	35	28	35	GY-07281-52		GY-07281-62	
50	45	34	45	GY-07281-54		GY-07281-64	
100	56	47	56	GY-07281-56		GY-07281-66	

Crucibles and Lids, Platinum

Don't oxidize with acids or salts

- Withstands temperatures to 2192°F (1700°C)

High temperature strength and high melting point, while prohibiting contamination at high temperatures up to 2192°F (1200°C).



Lids sold separately


Capacity (mL)	Top ID (mm)	Height (mm)	Crucibles		Covers	
			Catalog number	Price/ea†	Catalog number	Price/ea†
Standard-form crucibles						
10	25	28	GY-17830-15		GY-17831-15	
15	30	31	GY-17830-25		GY-17831-25	
20	33	35	GY-17830-35		GY-17831-35	
25	36	34	GY-17830-45		GY-17831-45	
30	40	36	GY-17830-55		GY-17831-55	
40	42	42	GY-17830-65		GY-17831-65	
50	45	44	GY-17830-75		GY-17831-75	
Wide-form crucibles						
50	53	29	GY-17830-82		GY-17831-82	
75	65	30	GY-17830-86		GY-17831-86	
100	70	33	GY-17830-88		GY-17831-88	
150	77	41	GY-17830-90		GY-17831-90	

†Platinum pricing may fluctuate with platinum market.

Crucible Tongs

Get a better grip

Nickel-plated steel tongs have positive-grip corrugated tips for safely handling dishes in muffle furnaces. Use with crucibles and dishes 1½" (3.8 cm) in diameter or wider.



Catalog number	Length	Price
GY-06443-50	20" (50.8 cm)	

Ceramic Marking Ink

Won't wear off

Use to permanently mark laboratory porcelain and ceramicware; apply with pen or brush. Markings become permanent when heated.



Catalog number	Bottle oz (mL)	Price
GY-17829-00	0.5 (15)	

Wait!

There's More at
ColeParmer.com

Find all the
product specs
you need!