

General-Purpose Latex Gloves

Better dexterity and elasticity

- Don't lose your grip—textured surface for better handling
- Get easy donning with beaded cuff and powdered version, or choose powder-free for reduced skin irritation
- Single-use, ideal for food and industrial applications
- Ambidextrous

NEW



Size	Length	Thickness	Qty/bx	Catalog number	Price per box
Powdered gloves					
Small	9½" (24 cm)	5.5 mil	100	GY-86232-35	
Medium				GY-86232-36	
Large				GY-86232-37	
X-large				GY-86232-38	
Powder-free gloves					
Small	9½" (24 cm)	5.5 mil	100	GY-86232-45	
Medium				GY-86232-46	
Large				GY-86232-47	
X-large				GY-86232-48	

Derma Thin® Latex Gloves

Get better comfort, fit, and grip

- Natural rubber formulation for better strength, durability, and wear
- Won't tear as easily
- Easy donning with lightly powered cornstarch coating
- Ambidextrous



Compliance: FDA 21 CFR parts 170–199.

Size	Length	Thickness	Qty/bx	Catalog number	Price per box
Small	9½" (24 cm)	5 mil	100	GY-86431-50	
Medium				GY-86431-51	
Large				GY-86431-52	

Powdered Latex Finger Cots

Provide protection while maintaining dexterity

Compliance: Meet ASTM D3772 for physical properties, manufactured under FDA standards.



Size	Length	Qty/bg	Catalog number	Price per bag
Small	2¾" (7 cm)	1440	GY-81605-32	
Medium			GY-81605-42	
Large			GY-81605-52	

Cole-Parmer® Nitrile Gloves

Don't lose your grip

- Soft feel with textured fingertips for better grip
- Powder-free
- Ambidextrous

NEW



Cole-Parmer®

Size	Length	Thickness	Qty	Catalog number	Price
Small	24.1 cm (9½")	4.0 mil	1 box of 100	GY-15954-42	
Medium				GY-15954-43	
Large				GY-15954-44	
X-large				GY-15954-45	
Small	24.1 cm (9½")	4.0 mil	10 boxes of 100 (1000 total)	GY-15954-24	
Medium				GY-15954-25	
Large				GY-15954-26	
X-large				GY-15954-46	

N-DEX® Nitrile Gloves

Keep your cool—special formulation for comfort, less fatigue

- Fewer punctures—three times the resistance of rubber or vinyl
- Ambidextrous

Puncture resistant gloves with-stand tough conditions



General and Laboratory Grade N-DEX Gloves protect your hands from the environment.

Industrial Grade N-DEX Plus Gloves are thicker for greater puncture and tear resistance. Gloves provide good resistance to chemicals and abrasion. Ideal for tough lab conditions, assembly, or manufacturing.

Compliance: FDA 21 CFR, parts 170-199, acceptable for food and pharmaceutical use.



CE

Size	Length	Thickness	Qty/pk	Catalog number	Price per pack
General and laboratory grade low-powder† N-DEX gloves					
Small	9½" (24 cm)	4 mil	100	GY-86287-70	
Medium				GY-86287-71	
Large				GY-86287-72	
X-large				GY-86287-73	
General and laboratory grade powder-free N-DEX gloves					
Small	9½" (24 cm)	4 mil	100	GY-86287-80	
Medium				GY-86287-81	
Large				GY-86287-82	
X-large				GY-86287-83	
Industrial grade N-DEX Plus low-powder† gloves					
Small	9½" (24 cm)	8 mil	50	GY-86287-60	
Medium				GY-86287-61	
Large				GY-86287-62	
X-large				GY-86287-63	

†USP grade cornstarch

Gloves

GREEN-DEX™ Biodegradable Nitrile Gloves**Maintain safety & the environment!**

- Reduce your carbon footprint—Eco Best Technology™ accelerates biodegradation
- Don't struggle—rolled cuff and smooth finish for easier donning
- Minimize contamination—disposable, single use
- Ambidextrous

**NEW**

Compliance: FDA 21 CFR for food contact and ASTM 5511 for biodegradation.



Size	Length	Thickness	Qty/bx	Catalog number	Price per box
X-small	9½" (24 cm)	4 mil	100	GY-86287-00	
Small				GY-86287-01	
Medium				GY-86287-02	
Large				GY-86287-03	
X-large				GY-86287-04	

Industrial-Grade Nitrile Gloves**Don't risk safety during tough jobs**

- Thicker for better protection
- Powder-free with beaded cuff for easier donning
- Powder-free
- Ambidextrous

**NEW**

Size	Length	Thickness	Qty/bx	Catalog number	Price per box
Small	9½" (24 cm)	5.5 mil	100	GY-86233-15	
Medium				GY-86233-16	
Large				GY-86233-17	
X-large				GY-86233-18	

Dura Flock™ Flock-Lined Nitrile Gloves**Keep hands dry**

- Moisture-absorbing lining keeps hands dry
- Superior protection with 8 mil thickness
- Powder-free textured gloves for better grip
- Ambidextrous



Size	Length	Thickness	Qty/bx	Catalog number	Price per box
Small	10½" (27 cm)	8 mil	50	GY-86484-40	
Medium				GY-86484-41	
Large				GY-86484-42	
X-large				GY-86484-43	

Touch N Tuff™ Nitrile Gloves**Keep dexterity without sacrificing safety**

- Thin formulation, yet resists punctures and tears
- No fillers—gloves are 100% nitrile
- Ambidextrous



Compliance: FDA regulations for food contact



Size	Length	Thickness	Qty/ bx	Catalog number	Price per box
Powdered gloves					
Small	9½" (24 cm)	4 mil	100	GY-81602-40	
Medium				GY-81602-50	
Large				GY-81602-60	
X-large				GY-81602-70	
Powder-free gloves					
Small	9½" (24 cm)	4 mil	100	GY-81602-45	
Medium				GY-81602-55	
Large				GY-81602-65	
X-large				GY-81602-75	

TNT™ Blue Disposable Nitrile Gloves**Durable and comfortable**

- Static dissipative to eliminate static electricity
- Textured fingertips for better gripping
- Powder-free
- Ambidextrous



Compliance: FDA regulations for food contact.

Size	Length	Thickness	Qty/bx	Catalog number	Price per box
Small	9½" (24 cm)	5 mil	100	GY-81601-04	
Medium				GY-81601-05	
Large				GY-81601-06	
X-large				GY-81601-07	

Glove Dispenser

Glove box dispenser holds one box of disposable gloves. Dispenser can be placed on benchtop or wall-mounted using hardware included. Measures (W x H x D) 10" x 5½" x 4" (25.2 x 14 x 10.2 cm).



[GY-63535-13](#) Glove box dispenser

Trilites® Gloves

Protect more area

- Longer, 10" gloves provide more protective area
- Better tear and chemical resistance
- Tri-polymer blend of neoprene, nitrile and rubber for more strength
- Powder-free
- Ambidextrous



Compliance: USDA accepted per FDA 21 CFR, parts 170-199.

Size	Length	Thickness	Qty/box	Catalog number	Price per box
Small	10" (25 cm)	6 mil	100	GY-40215-00	
Medium				GY-40215-01	
Large				GY-40215-02	
X-large				GY-40215-03	

Nitrile Finger Cots

Reduce finger fatigue

- Nitrile formulation offers superior stretch and comfort
- Latex-free to eliminate reactions
- Antistatic—eliminates damaging static discharge
- Powder-free



Size	Length	Thickness	Qty/bag	Catalog number	Price per bag
Medium	2 3/4"	3 mil	720	GY-81605-92	
Large	(7 cm)			GY-81605-94	

Chemical-Resistant PVA Gloves

Will not swell or shrink

- Reduced absorption to eliminate swelling and shrinking
- Provide extreme resistance to strong solvents
- More comfortable than neoprene-coated gauntlets, weighing half as much
- Resist snags, punctures, abrasions and cuts



Size	Length	Qty/pk	Catalog number	Price per pack
9	14" (35 cm)	12 pair	GY-09715-04	
10			GY-09715-05	

Sol-Vex® Nitrile Gloves

Don't risk chemical burns



- High-performance nitrile compound provides outstanding strength and chemical resistance
- Choose lined or unlined versions

Compliance: Made with FDA compliant materials



Size	Length	Thickness	Qty/pk	Catalog number	Price per pack
Unlined gloves					
9	13"	11 mil	12 pairs	GY-09715-30	
10	(33 cm)			GY-09715-40	
9	18"	22 mil		GY-81603-40	
10	(46 cm)			GY-81603-60	
Lined gloves					
9	13"	15 mil	12 pairs	GY-09715-32	
10	(33 cm)			GY-09715-42	

DON'T RISK IT!

Make sure your gloves are compatible with your materials.

ColeParmer.com/SafetyChemGuide



Silver Shield® Gloves and Sleeves

Protect against aggressive chemicals



Gloves 86358-02 shown with sleeves 86358-03

- Resist permeation by a wide range of aromatic, chlorinated, aliphatic solvents, ketones, and esters
- Ideal for lab or industrial use
- Ambidextrous and available in three sizes
- Sleeves feature an elastic cuff at both ends

Description	Size	Qty/pk	Catalog number	Price per pack
Gloves	Small	10 pairs	GY-86358-00	
	Medium		GY-86358-01	
	Large		GY-86358-02	
Sleeves	18" (46 cm) L	1 pair	GY-86358-03	

Gloves

Butyl II Gloves**Stay safe during clean-ups**

- Non-bulky and non-restrictive for easier handling
- Excellent chemical resistance
- Long, 14" gloves for added protection
- Ideal for use with hazmat suits

Compliance: ASTM F739.



Size	Description	Catalog number	Price per pair
Medium	Smooth gloves, 14-mil thick	GY-86357-01	
Large		GY-86357-02	
X-large		GY-86357-04	
Medium	Textured gloves, 14-mil thick	GY-86357-11	
Large		GY-86357-12	
X-large		GY-86357-13	
Medium	Heavyweight smooth gloves, 30-mil thick	GY-86357-20	
Large		GY-86357-21	
X-large		GY-86357-22	

Viton® Gloves**Get flexible, heavy-duty protection**

- Superior resistance to PCBs, as well as chlorinated and aromatic solvents



Size	Length	Thickness	Catalog number	Price per pair
Small	14" (35 cm)	12 mil	GY-86288-20	
Medium			GY-86288-21	
Large			GY-86288-22	
X-large	14" (35 cm)	30 mil	GY-86288-23	
Large			GY-86288-32	
X-large			GY-86288-33	

Cleanroom Gloves**Eliminate contamination**

- Low residues including particulates, ionic extractables, and nonvolatiles
- Textured surface for better gripping
- Latex-free to eliminate adverse reactions
- Surface resistivity of 10¹¹ ohm/square for static decay
- Ambidextrous

Compliance:
ISO Class 4 or 5



Size	Length	Thickness	Qty/bg	Catalog number	Price per bag
Latex ISO Class 4 cleanroom gloves					
Small	12" (30 cm)	6.3 mil	100	GY-86313-35	
Medium				GY-86313-36	
Large				GY-86313-37	
X-large				GY-86313-38	
Nitrile ISO Class 4 cleanroom gloves					
Small	12" (30 cm)	5.2 mil	100	GY-86313-45	
Medium				GY-86313-46	
Large				GY-86313-47	
X-large				GY-86313-48	
Latex ISO Class 5 cleanroom gloves					
Small	12" (30 cm)	6.3 mil	100	GY-86313-15	
Medium				GY-86313-16	
Large				GY-86313-17	
X-large				GY-86313-18	
Nitrile ISO Class 5 cleanroom gloves					
Small	12" (30 cm)	5.2 mil	100	GY-86313-25	
Medium				GY-86313-26	
Large				GY-86313-27	
X-large				GY-86313-28	

Vinyl Cleanroom Gloves**Don't contaminate samples**

- Latex-free, PVC copolymer is naturally antistatic
- Micro-textured surface for better gripping
- Ambidextrous
- Packed 100 gloves per polybag

Compliance:
ISO Class 5, USDA,
and IEST-RP-CC-005.2



Size	Length	Thickness	Qty/bg	Catalog number	Price per bag
Medium	12" (30 cm)	6 mil	100	GY-86580-21	
Large				GY-86580-22	
X-large				GY-86580-23	

Biohazard Autoclave Gloves

Don't get burned

- Heat resistant to 450°F (232°C)—ideal for removing hot objects from an oven
- Machine washable
- One size fits all



Length	Catalog number	Price per pair
13" (33 cm)	GY-06496-40	
18½" (47 cm)	GY-06496-42	

Cotton Inspection Gloves

Stop leaving fingerprints

- Minimize cross-contamination from human hands
- Washable and reversible for reusing
- Ideal for liner between liquid-proof gloves and hands
- Ambidextrous



Size	Qty/pk	Catalog number	Price per pack
Lightweight gloves			
Women's	12 pairs	GY-86369-00	
Men's		GY-86369-02	
Medium-weight gloves			
Men's	12 pairs	GY-86369-05	

HyFlex® Lite Precision Work Gloves

Easily handle small parts

- Stretch nylon liner for better tactile sensitivity



Size	Catalog number	Price per pair
7	GY-09719-00	
8	GY-09719-01	
9	GY-09719-02	
10	GY-09719-03	



HyFlex™ Work Gloves

Maintain your grip

- Nylon knit gloves provide a snug yet cool fit for maximum dexterity and tactile sensitivity
- Ideal for general assembly and handling small parts
- Nitrile foam coating resists snags and abrasions
- Excellent grip in light oil applications



Size	Catalog number	Price per pair
7	GY-81611-30	
8	GY-81611-32	
9	GY-81611-34	

Nitri-Knit® Extended Length Gloves

Protect your arms

- Extended length provides protection of upper arms
- Won't peel apart like glued-in liner gloves—nitrile is permanently bonded to inner liner
- Won't swell, sag, or weaken from fat absorption
- Textured finish for better gripping
- Hand portion is treated with antibacterial and antimildew agent
- Launderable
- Latex-free



Size	Length	Catalog number	Price per pair
7	24" (61 cm)	GY-86911-00	
8		GY-86911-01	
9		GY-86911-02	
10		GY-86911-03	

Glove Sizing Guide

Measure the circumference of your hand around the palm area with fingers together and hand relaxed. This measurement in inches is closest to your actual glove size.

For example, 7" is equal to a size 7 glove. Sizes may vary between different manufacturers and styles. Except where mentioned, gloves are sized according to men's hands. For easy size comparison, see table at right.



Glove Size Comparison

Numerical size	6-7	7-8	8-9	9-10	11
General size	XS	S	M	L	XL

Teky's Tips



Gloves

Thermaprene™ Heat-Resistant, Lined Neoprene Gloves

Don't give up comfort for safety

- Soft thermal liner provides heat resistance and comfort
- Chemical-resistant neoprene coating
- Protects against punctures
- Nonslip grip



Size	Length	Catalog number	Price per pair
8	26" (66 cm)	GY-81609-04	
10		GY-81609-06	

Temp-Tec® Heat-Insulated Neoprene Gloves

Don't wonder when to replace

- Two-tone indicator alerts you when to replace gloves
- Good chemical resistance
- Roughened surface for better gripping
- Withstand temperatures from –100 to 480°F (–73 to 249°C)

Compliance: FDA 21 CFR 177.2600.



Insulated neoprene gloves give you excellent protection from extreme hot or cold.



Size	Length	Thickness	Catalog number	Price per pair
Medium	14" (35 cm)	Heavy gauge	GY-40211-30	
Large			GY-40211-31	
X-large			GY-40211-32	

Temp-Proof® 333 Heat-Insulated Nitrile Gloves

Protect from both hot and cold

- Thermal insulation from –22 to 480°F (–30 to 249°C)
- Launderable to reduce costs and environmental waste



Size	Length	Thickness	Catalog number	Price per pair
9	12" (30 cm)	Heavyweight	GY-40211-80	

Extreme-Temperature Gloves
Stand up to the heat!

- Withstand dry temperatures from –210 to 1400°F (–134 to 760°C)
- Won't char or leave residue
- Free of PCB, asbestos, and fiberglass
- Will not leave particulates behind—flexible nylon lining



Compliance: ISO Class 5; ESD per DODD-HDBCV-263 and MIL-B-81705.

Size	Length	Catalog number	Price per pair
Max. temperature 1000°F (537°C) gloves			
Large	14" (35 cm)	GY-33689-20	
Max. temperature 1400°F (760°C) gloves			
Large	14" (35 cm)	GY-33689-21	
	18" (46 cm)	GY-33689-22	
	27" (69 cm)	GY-33689-23	

High-Temperature Gloves
Won't let you get burned

- Material resistant up to 2000°F (1093°C)
- Asbestos-free, cotton-backed wool liner
- Extra strong with aramid thread throughout
- One size fits all
- Made in USA

Compliance: European EN407 industrial standard.



86364-02



86364-10



Description	Peak temperature†	Catalog number	Qty	Price
Zetex glove, 14"L (35 cm)	1000°F (538°C)	GY-86364-00	1 pair	
ZetexPlus glove, 14"L (35 cm)	2000°F (1093°C)	GY-86364-10	1 pair	
Zetex finger cots	1000°F (538°C)	GY-86364-02	12/pk	
Aluminized Zetex glove, 14"L (35 cm)	1000°F (538°C) radiant heat	GY-86364-20	1 pair	

†Peak temperature is dependent on duration, application, and environment.

Cryogenic Gloves

Keep hands from freezing

- Protect as low as –256°F (–160°C)
- Maximum protection while remaining flexible

Water-Resistant Gloves provide protection in closed systems. A stay-dry lining keeps you warm and comfortable. **Waterproof Gloves** have a nylon exterior and a seamless waterproof lining. Gloves keep you dry and wick moisture away from your hands. **Cryo-Industrial® Gloves** feature a rugged, abrasion-resistant outer shell that protects the glove and provides an excellent grip for heavy-duty applications. Waterproof liner protects against cryogenic splashes while the inner lining wicks moisture from your hands for enhanced comfort.

Note: Do not immerse gloves in liquid nitrogen.



Description	Size	Length	Catalog number	Price per pair
Water-resistant Cryo-Glove gloves				
Wrist	Medium	12" (30 cm)	GY-09113-02	
	Large		GY-09113-04	
Mid-arm	Medium	15" (38 cm)	GY-09113-12	
	Large		GY-09113-14	
	X-large		GY-09113-16	
Elbow	Medium	19" (48 cm)	GY-09113-22	
	Large		GY-09113-24	
	X-large		GY-09113-26	
Waterproof Cryo-Glove gloves				
Mid-arm	Medium	15" (38 cm)	GY-09113-30	
	Large		GY-09113-32	
	X-large		GY-09113-34	
Elbow	Medium	19" (48 cm)	GY-09113-35	
	Large		GY-09113-37	
	X-large		GY-09113-39	
Cryo-Industrial heavy-duty gloves				
Mid-arm	Medium	15" (38 cm)	GY-09113-62	
	Large		GY-09113-64	
	X-large		GY-09113-66	

Drybox Gloves

Protect your entire arm

- Shoulder-length, seamless protection

Neoprene Gloves feature high tensile strength and flexibility for work that requires high sensitivity and good dexterity.

Hypalon® Gloves offer outstanding abrasion protection and resistance to ozone or oxidizing chemicals. Natural softness improves worker comfort.

Butyl Gloves allow highest impermeability to water, vapors and gases. They maintain flexibility even at low temperatures and are resistant to oxygenated solvents and most oxidizing chemicals.



Material	Size	Length	Thickness	Catalog number	Price per pair
Neoprene	9¾	32" (81 cm)	15 mil	GY-86569-00	
			30 mil	GY-86569-02	
Hypalon	9¾	32" (81 cm)	15 mil	GY-86569-10	
			30 mil	GY-86569-12	
Butyl	9¾	32" (81 cm)	15 mil	GY-86569-20	
			30 mil	GY-86569-22	

Insulating Rubber Glove Kits

Don't compromise safety when working with electrical power

- Includes everything you need—insulating rubber gloves, appropriate leather protector, and storage bag
- Protect gloves from cuts, abrasions, and punctures with included leather protector

Tip: In accordance with NFPA 70E, anything over 50 volts requires insulating rubber gloves.



Kits include leather protector, storage bag and insulating rubber gloves



Class	Length	Max use voltage AC / DC	Size	Red gloves		Black gloves		Yellow gloves	
				Catalog number	Price	Catalog number	Price	Catalog number	Price
00	11" (279 mm)	500 / 750	8	GY-86435-00		GY-86435-05		—	—
00	11" (279 mm)	500 / 750	9	GY-86435-01		GY-86435-06		—	—
00	11" (279 mm)	500 / 750	10	GY-86435-02		GY-86435-07		—	—
00	11" (279 mm)	500 / 750	11	GY-86435-03		GY-86435-08		—	—
00	11" (279 mm)	500 / 750	12	GY-86435-04		GY-86435-09		—	—
0	11" (279 mm)	1000 / 1500	8	GY-86435-10		GY-86435-15		GY-86435-20	
0	11" (279 mm)	1000 / 1500	9	GY-86435-11		GY-86435-16		GY-86435-21	
0	11" (279 mm)	1000 / 1500	10	GY-86435-12		GY-86435-17		GY-86435-22	
0	11" (279 mm)	1000 / 1500	11	GY-86435-13		GY-86435-18		GY-86435-23	
0	11" (279 mm)	1000 / 1500	12	GY-86435-14		GY-86435-19		GY-86435-24	
2	14" (356 mm)	17,000 / 25,500	8	—	—	GY-86435-25		—	—
2	14" (356 mm)	17,000 / 25,500	9	—	—	GY-86435-26		—	—
2	14" (356 mm)	17,000 / 25,500	10	—	—	GY-86435-27		—	—
2	14" (356 mm)	17,000 / 25,500	11	—	—	GY-86435-28		—	—
2	14" (356 mm)	17,000 / 25,500	12	—	—	GY-86435-29		—	—