

Ice

Flaked-Ice Maker/Dispenser

Eliminate cross-contamination

- Touch-free actuation provides convenient one-handed operation
- Up to 25% less fan noise
- Ice only or ice/water selectable
- Easy access for most containers



Condenser unit	Ice/24 hr	Bin capacity	Power (VAC, Hz)	Catalog number	Price
Air cooled	392 lb (178 kg)	12 lb (5.4 kg)	115, 60	GY-55100-06	

Wait!

There's More at ColeParmer.com

Find more than **100,000** products to choose from!

Junior Dry Ice Maker

No more wasted dry ice

- Produces only what you need in brick or snow form
- Simply attach to compressed LCO₂ tank, fill bag, and squeeze if hard pack desired
- Makes -0.77 kg of flakes per bag

Note: LCO₂ Cylinder must have a dip tube.

What's included: collection bag, clamp, and plastic cone.



Description	Catalog number	Price
Junior dry ice maker	GY-03716-32	

[GY-03716-90](#) Extra bags. Pack of 3

Cubed- and Flaked-Ice Machines

Don't run out of ice

- High-volume production up to 395 lbs per day
- Convenient size fits under standard lab benches
- Get easy access through slide-away door

Cubed-Ice Makers are constructed of durable rust-free rotocast plastic for greater insulation and reduced melting of ice.

Flaked-Ice Makers produce up to 395 lbs of flaked ice per day and store 140 lbs in the self-contained bin. Flaked ice quickly conforms to your container or sample tubes for greater cooling capabilities.



Cubed ice



Flaked ice



Condenser unit	Ice/24 hr	Bin capacity	Power (VAC, Hz)	Catalog number	Price
Cubed-ice makers					
Air cooled	84 lb (38 kg)	24 lb (11 kg)	115, 60	GY-55100-05	
Air cooled	150 lb (68 kg)	86 lb (39 kg)	115, 60	GY-55100-01	
Water cooled				GY-55100-02	
Flaked-ice makers					
Air cooled	395 lb (179 kg)	80 lb (36 kg)	115, 60	GY-55100-03	
Water cooled				GY-55100-04	

Dry Ice Machine

Get more dry ice out of your CO₂

- Provides 30% more yield from standard 50-lb CO₂ cylinder
- No need to dismantle to get dry ice
- No tip, stainless steel design
- Makes compact 1-lb blocks that measure 6"L x 4"W x 1 1/2"H

What's included: CO₂ supply regulator valve, hose, and fitted brass safety release valve that automatically vents excess gas pressure.



Description	Catalog number	Price
Dry ice machine	GY-03716-03	
Dry ice machine for European cylinders	GY-03716-04	

Laptop Coolers

Say goodbye to messy ice baths

- Body is molded polycarbonate filled with a nontoxic insulating gel
- Clear polycarbonate gel-filled white polycarbonate lid on some models
- Lids feature a write-on surface for sample identification

Models 44001-00, -10, and -25 include adapters to accommodate 0.5-mL tubes. Models 44001-20 and -30 hold tubes up to 125-mm tall when lid is latched. All models have a stainless steel handle for a secure hold during transportation. Stackable design saves space.

Tube size	Tubes held	Temperature maintained	Dimensions (L x W x H)		Catalog number	Price
			in.	cm		
0.5 to 2.0 mL	12	–15°C, 1 hour	5 ⁵ / ₁₆ x 4 ¹ / ₄ x 4 ¹⁵ / ₁₆	13.5 x 10.8 x 12.5	GY-44001-00	
0.5 to 2.0 mL	32	–15°C, 2 hour	9 ⁹ / ₁₆ x 6 ³ / ₁₆ x 5 ³ / ₄	24.3 x 15.7 x 14.6	GY-44001-10†	
0.5 to 2.0 mL	32	1°C, 5 hour	9 ⁹ / ₁₆ x 6 ³ / ₁₆ x 5 ³ / ₄	24.3 x 15.7 x 14.6	GY-44001-25†	
12 to 13 mm dia	12	1°C, 5 hour	7 ³ / ₄ x 5 ¹ / ₂ x 7 ¹ / ₂	19.7 x 14.0 x 19.1	GY-44001-20	
16 to 17 mm dia	12	1°C, 5 hour	7 ³ / ₄ x 5 ¹ / ₂ x 7 ¹ / ₂	19.7 x 14.0 x 19.1	GY-44001-30	

†Coolers have a gel-filled white polycarbonate lid.



Polyurethane Foam Ice Bucket with Lid

Don't let ice melt

- Store up to 4 liters of ice for extended periods of time
- Insulation prevents ice from melting for hours
- Tapered design for convenient stacking
- Waterproof—will not rust, leak, or sweat



Bucket is inert to acetone, liquid nitrogen, dry ice, and salt.

Color	Capacity	Dimensions (H x dia)		Catalog number	Price
		in.	cm		
Black	4 L	7 ³ / ₄ x 9	18.7 x 22.9	GY-06274-21	

Portable Coolers

Keep enzymes or samples cold for transport

- Refrigerate overnight to maintain specified temperatures
- Small enough to fit in your hand



Tubes held x tube size	Temperature maintained	Catalog number	Price
12 x 1.5-mL or 12 x 2-mL	0°C, 10 hrs	GY-44540-00	
	15°C, 5 hrs	GY-44540-02	

Freezers

Keep temperatures constant

- Designed for routine lab use
- Manual defrost

Undercounter Freezer features an easy-to-clean plastic interior, two fixed wire shelves, and three door shelves. Units are UL listed.

Upright Freezer is ideal for lab applications where storage of noncritical materials is important. Easy-to-clean painted interior features four fixed position wire shelves and six molded door shelves. Keylocked, hinged door provides added security.



44202-10

Capacity	Temperature range	Power		Catalog number	Price
		VAC	Hz		
Undercounter freezer					
5.6 cu ft (142 L)	–20 to –15°C (–4 to 5°F)	120	60	GY-44202-10	
Upright freezer					
20 cu ft (566 L)	–20 to –12°C (–4 to 10°F)	120	60	GY-44202-40	
		240	50	GY-44202-45	

Wait!

There's More at
ColeParmer.com

Find all the
product specs
you need!