


Absorbent Lab Soakers

Protect glassware and surfaces

- No more difficult clean-ups with quick absorbing material
- Skid resistant, polyethylene backing grips surface and retains liquids



Sq ft/pk	Size (W x L)	Catalog number	Qty/pk	Price/pk
Absorbent lab soakers, absorb 750 mL/m²				
1000	20" x 300 ft (50.8 cm x 91.4 m) roll	GY-06803-30	2 rolls	
875	18" x 20" (45.7 x 50.8 cm) mat	GY-06803-32	350 mats	
Super absorbent lab soakers, absorb 1050 mL/m²				
417	20" x 250 ft (50.8 cm x 76.2 m) roll	GY-06803-40	1 roll	
125	18" x 40" (45.7 x 101.6 cm) mat	GY-06803-42	25 mats	
150	20" x 43" (50.8 x 109.2 cm) mat	GY-06803-44	25 mats	

Surface Protectors

Won't let corrosive materials through

- Protect surfaces from corrosion and absorption of toxic materials
- Pressure-sensitive backing—just peel backing and apply to surface
- Just wipe clean or wash down




FEP-coated is ideal for use in kitchens and laboratories.

Support material	Support material color	Thickness	Max temp.	Roll width	Roll length	Catalog number	Price/roll
Vinyl	White	8 mm	200°F (93°C)	12½" (31.8 cm)	15 ft (4.5 m)	GY-06804-20	
				25" (63.5 cm)		GY-06804-00	
Aluminum	Silver	2 mm	360°F (182°C)	12½" (31.8 cm)		GY-06804-30	
				25" (63.5 cm)		GY-06804-10	

Liners and Sheets

Clean up spills easily

- Disposable liners absorb most liquids
- Easily cut to fit any area
- Won't leak or skid with polyethylene backing



Color	Size (W x L)	Catalog number	Qty	Price
White	20" x 50 ft (50.8 cm x 15 m)	GY-06805-01	1 roll	
Orange	20" x 50 ft (50.8 cm x 15 m)	GY-06805-03	1 roll	
White	20" x 24" (50.8 cm x 61.0 cm)	GY-06805-05	100 sheets/case	

Absorbent Bench Protectors

No need for scissors

- Precut sheets for quicker placement
- Keep spills contained with absorbent upper surface and polycoated base
- Absorb twice as much with BenchGuard Extra



Description	Absorbency	Size (W x L)	Catalog number	Pack	Price
BenchGuard	400 mL/m ²	19.3" x 164 ft (49 cm x 50 m)	GY-06805-00	1 roll	
		19.3" x 24" (49 cm x 61 cm)	GY-06806-00	100 sheets	
BenchGuard Extra	800 mL/m ²	19.3" x 164 ft (49 cm x 50 m)	GY-06805-10	1 roll	
		19.3" x 24" (49 cm x 61 cm)	GY-06806-10	50 sheets	

Liners, Neoprene-Coated Fiber

Contour to your needs

- Cut to fit—easily contours to irregular surfaces
- Deadens sounds of items being placed or dropped
- Autoclavable



(Shown at actual size)

Mesh opening	Material	Thickness	Max temp.	Roll size (W x L)	Catalog number	Price/roll
¼" (6 mm)	Black	0.125"	250°F (121°C)	26" x 50 ft (66.0 cm x 15 m)	GY-09400-01	
	Red				GY-09400-02	
	Green				GY-09400-03	
	Almond				GY-09400-04	

Protective Netting, HDPE

Minimize impact on glassware

- Simply place over glassware to protect from bumps or flying debris
- Cut to fit your needs

Approximate diameter	Color	Thickness	Maximum temperature	Roll length	Catalog number	Price/roll
Close-mesh netting						
3/8" to 5/8" (1.0 to 1.6 cm)	Orange	0.050"	175°F (80°C)	820 ft (250 m)	GY-09405-15	
5/8" to 1" (1.6 to 2.5 cm)	Blue	0.050"			GY-09405-21	
1" to 2" (2.5 to 5.1 cm)	Yellow	0.050"			GY-09405-25	
2" to 4" (5.1 to 10.2 cm)	Red	0.070"	175°F (80°C)	164 ft (50 m)	GY-09405-30	
4" to 6" (10.2 to 15.2 cm)	Green	0.070"			GY-09405-35	
6" to 8" (15.2 to 20.3 cm)	Brown	0.070"	175°F (80°C)	164 ft (50 m)	GY-09405-40	
8" to 10" (20.3 to 25.4 cm)	Crystal	0.070"			GY-09405-45	
Heavy-gauge netting —more rigid for extra protection						
1 1/4" to 1 5/8" (3.2 to 4.1 cm)	Yellow	0.110"	175°F (80°C)	82 ft (25 m)	GY-09405-50	
2 1/2" to 3" (6.4 to 7.6 cm)	Green	0.110"			GY-09405-55	



Close mesh



Heavy-gauge

(Shown at actual size)

Lab Matting, PTFE Fluoropolymer

Don't let items go down the drain

- Line sinks and shelves to minimize impact
- Easily cut to size
- 1/4" (6 mm) mesh opening
- Autoclavable



(Shown at actual size)

Material	Thickness	Maximum temperature	Roll size (W x L)	Catalog number	Price/roll
White	0.0625"	500°F (260°C)	12" x 10 ft (31 cm x 3 m)	GY-09406-00	

Cole-Parmer® PTFE Film

Stain and stick resistance

- Protect surfaces from chemicals and high temperatures
- Simply wipe away spills
- Cut to size with scissors for custom fit
- Reusable




Cole-Parmer®

Color	Thickness	Maximum temperature	Roll size (W x L)	Catalog number	Price/roll
Opaque	0.005" (1.3 mm)	500°F (260°C)	12" x 25 ft (30.5 cm x 7.6 m)	GY-08277-15	

Ribbed Matting, LDPE

Minimize broken glass

- Prevent glass breakage by reducing impact on surfaces
- Solvent resistant



Color	Thickness	Maximum temperature	Roll size (W x L)	Catalog number	Price/roll
Clear	0.036"	176°F (80°C)	30" x 15 ft (76.2 cm x 4.6 m)	GY-06809-10	
			30" x 50 ft (76.2 cm x 15 m)	GY-06809-20	

Don't Forget...

Cleanroom Wipes and Mats
See pages 204–210

