

Scoops & Spoons

Round Scoops, Aluminum

Handle won't separate from scoop—durable one-piece construction

- Save money with aluminum compared to stainless steel
- Don't lose your grip—hold tight with contoured handle



Capacity mL (oz)	Scoop (L x W)	Catalog number	Qty/ pk	Price/ pk
150 (5)	4 3/4" x 2 3/4" (12.1 x 6.0 cm)	GY-07207-00	1	
350 (12)	6 1/4" x 3 1/4" (15.9 x 8.3 cm)	GY-07207-02	1	
700 (24)	8" x 4" (20.3 x 10.2 cm)	GY-07207-04	1	
1100 (38)	9" x 4 3/4" (22.9 x 12.1 cm)	GY-07207-06	1	
1700 (57)	10 1/2" x 5 1/2" (26.7 x 14.0 cm)	GY-07207-08	1	
2500 (84)	12" x 6 1/4" (30.5 x 15.9 cm)	GY-07207-12	1	

Cole-Parmer® Scoops, SS

Won't corrode—ideal for food-grade sampling

- Choose 304 stainless steel for higher corrosion resistance



Cole-Parmer®

Capacity mL (oz)	Scoop (L x dia)	Catalog number	Qty/ pk	Price/ pk
201 Stainless steel scoops				
175 (5.5)	5 1/2" x 3" (14 x 7.6 cm)	GY-07205-01	1	
1550 (52)	10" x 5 3/4" (25.4 x 14.6 cm)	GY-07205-03†	1	
304 Stainless steel scoop				
1550 (52)	7 1/8" x 5 1/4" (18.1 x 13.3 cm)	GY-07210-00	1	

†Model 07205-03 has a hook on the end for convenient storage.

Square Scoop, PE

Stands on end—will not tip over

- Get large quantity scooping easily



Capacity mL (oz)	Dimensions (L x W x H)	Catalog number	Qty/ pk	Price/ pk
2000 (64)	8 1/4" x 6" x 5" (21.0 x 15.2 x 12.7 cm)	GY-06277-00	1	

Round Scoops, HDPE

Scoop easier with flat bottom

- Get a better grip with 4 1/2" (11.4 cm) handle



Capacity mL (oz)	Scoop (L x W)	Catalog number	Qty/ pk	Price/ pk
350 (12)	8" x 4 1/4" (20.3 x 10.8 cm)	GY-06282-20	1	
500 (16)	8" x 6" (20.3 x 15.2 cm)	GY-06282-30	1	
1000 (32)	8" x 6" (20.3 x 15.2 cm)	GY-06282-40	1	

Rectangular Scoops, PC

Withstands impact

- Made of tough, transparent polycarbonate
- Store out of the way with convenient handle hole for hanging



Capacity mL (oz)	Scoop (L x W)	Catalog number	Qty/ pk	Price/ pk
200 (6)	5" x 3" (12.7 x 7.6 cm)	GY-66600-00†	1	
350 (12)	6" x 3" (15.2 x 7.6 cm)	GY-66600-10	1	
1000 (32)	7 1/4" x 4 3/4" (18.4 x 12.1 cm)	GY-66600-20†	1	
2000 (64)	9 1/4" x 5 7/8" (23.5 x 14.9 cm)	GY-66600-30†	1	

†These items are USDA-accepted and comply with FDA regulations.

Sterileware® Sense-able™ Detectable Scoops, PS

Prevent product recalls

- Made of metal and X-ray equipment detectable material
- FDA-grade resin is detectable in small fragments
- Individually gamma sterilized and packaged



Capacity mL (oz)	Scoop (L x W x H)	Catalog number	Qty/ bx	Price/ bx
60 (2)	3" x 1 1/2" x 1" (7.6 x 3.8 x 2.5 cm)	GY-06263-65	100	
125 (4)	3 1/2" x 1 3/4" x 1 1/4" (8.9 x 4.4 x 3.2 cm)	GY-06263-70	100	
250 (8)	5" x 2" x 2" (12.7 x 5.1 x 5.1 cm)	GY-06263-75	100	

Sterile Sampling Scoops, PS

No need to presterilize

- Come sterilized in individual polyethylene bags
- Ideal for single use



Capacity mL (oz)	Scoop (L x W)		Handle length		Catalog number	Qty/ cs	Price/ cs
	in.	cm	in.	cm			
60 (2)	3 x 1	7.6 x 2.5	4	10.2	GY-06260-10	100	
120 (4)	3½ x 1	8.9 x 2.5			GY-06260-20		

Sterile Sampling Spatula, PS

Get more volume per scoop

- Designed to hold twice the volume of standard spatulas
- Break up caked material easily with sturdy edge
- Individually sterilized and bagged



Overall length	Catalog number	Qty/cs	Price/cs
9" (22.9 cm)	GY-06262-00	100	

Cole-Parmer® Double-Ended Spoon Set, 304 SS

Sample from both ends

- Minimize discomfort—solid rounded shafts for more comfortable handling

Set contains three spoons: a 5⅞" (14.9 cm) long spoon with 0.5 mL and 2.0 mL bowls, an 8¼" (21 cm) long spoon with 1.6 mL and 5.0 mL bowls, and an 11⅞" (30.2 cm) long spoon with 2.0 mL and 8.0 mL bowls.



Set contains three spoon sizes



Qty/set	Catalog number	Price/set
Three spoons	GY-06261-00	

Autoclavable Scoops, PP

Eliminate contamination risks

- Autoclavable at 250°F (121°C) to maintain sample integrity
- Rigid, strong design stands up to your needs



Capacity mL (oz)	Scoop (L x W)		Catalog number	Qty/ pk	Price/ pk
	in.	cm			
22 (¾)	2⅜ x 1⅛	6.0 x 2.9	GY-06275-15	12	
60 (2)	3 x 1½	7.6 x 3.8	GY-06275-16	12	
75 (2½)	3½ x 1½	8.9 x 3.8	GY-06275-17	12	
100 (3½)	3⅝ x 1⅞	9.2 x 4.8	GY-06275-18	12	
200 (6¾)	5 x 2	12.7 x 5.1	GY-06275-19	12	
1100 (37)	8 x 4¾	20.3 x 12.1	GY-06275-20	6	

Sterile Sampling Spoons, PS

Stop washing spoons

- Individually sterilized and bagged in a polyethylene bag
- Ideal for single use




Capacity	Handle length	Catalog number	Qty/cs	Price/cs
¼ tsp (1.2 mL)	6⅞" (15.5 cm)	GY-06263-10	200	
½ tsp (2.5 mL)	5⅞" (14.9 cm)	GY-06263-20	200	
1 tsp (4.9 mL)	5½" (14.0 cm)	GY-06263-30	200	
1 tbsp (14.8 mL)	5¼" (13.3 cm)	GY-06263-40	200	

Lab Spoon/Spatula, 304 SS

Easier sampling from large containers

- Total length of 12" (31 cm)



Spoon capacity	Spoon end	Spatula end	Catalog number	Qty/pk	Price/pk
5 mL (0.18 oz)	1½" (3.8 cm) W	1⅜" (3 cm) W	GY-01019-14	1	