

Single Test Verification Materials (STVM)

- Gas oil with a certified value provided for a single test method
- Characterised using a network of competent laboratories to provide international traceability
- Use this standard to cross-check instrument performance and to assist in operator training with a typical material
- Small pack size to avoid cross-contamination and help meet verification requirements for audited laboratories and those which comply with ISO 17025
- 2 year shelf life from manufacture

Single Test Verification Material - Gas Oil 50 mL

- 50 mL
- Produced in general accordance with the principles of ISO 17034 and ISO 33405 with documentation compliant to ISO 33401



Part Number:	Test Name:	Test Method:	Range:
99867-0	Sulfur WD X-ray	ASTM D2622; ISO 20884	0.001 to 0.013%(m/m)
99869-0	Polycyclic Aromatics HPLC	ASTM D6591-IP 391; EN 12916	0.6 to 3.6%(m/m)
99981-0	Sulphur	ASTM D5453; IP 490; ISO 20846	3 to 10 mg/kg

SIMDIS Single Test Verification Material - Gas Oil 10 mL (99906-0)

- 10 mL
- Produced in general accordance with the principles of ISO 17034 and ISO 33405 with documentation compliant to ISO 33401



Test Name:	Test Method:	Range:	Amount/Test:
SIMDIS IBP	ASTM D2887-IP 406	120 to 190 °C	1 mL
SIMDIS 10%	ASTM D2887-IP 406	177 to 242 °C	1 mL
SIMDIS 50%	ASTM D2887-IP 406	260 to 290 °C	1 mL
SIMDIS 90%	ASTM D2887-IP 406	320 to 355 °C	1 mL
SIMDIS 95%	ASTM D2887-IP 406	335 to 375 °C	1 mL
SIMDIS FBP	ASTM D2887-IP 406	350 to 414 °C	1 mL